Table of Contents:

English (en)	page 3
Czech (cs)	page 70
Danish (da)	page 139
German (de)	page 207
Greek (el)	page 281
Spanish (es)	page 354
Finnish (fi)	page 424
French (fr)	page 495
Icelandic (is)	page 566
Italian (it)	page 632
Dutch (nl)	page 702
Norwegian (no)	page 770
Polish (pl)	page 838
Portuguese (pt)	page 905
Swedish (sv)	page 974
Turkish (tr)	page 1043

Summary of Safety and Clinical Performance (SSCP)

Arthrosamid®

Document ID		018.SSCP.00023.03	
Version no.	001	Document date	2024-06-04

Table of contents

Abbreviations	5
Summary of safety and clinical performance - Part I - Professional users	6
1. Device identification and general information	6
1.1 Device trade name	6
1.2 Manufacturer	6
1.3 Manufacturer single registration number (SRN)	6
1.4 Basic UDI-DI	6
1.5 Medical device nomenclature description	6
1.6 Class of device	7
1.7 Year when the first CE certificate was issued covering the device	7
1.8 Notified body name and identification number	7
2. Intended purpose, indication, target population and contraindications	7
2.1 Intended purpose, indication and target population	7
2.2 Contraindications and/or limitations	7
3. Device description	8
3.1 Description of device	8
3.1.1 Mode of action	8
3.2 Previous generations of the device	8
3.3 Any other devices or products intended to be used in combination with the device	8
4. Residual risks and warnings	8
4.1 Residual risks and undesirable effects	8
4.2 Warnings and precautions	9
4.2.1 Warnings	9
4.2.2 Precautions	9
4.3 Other information	10
5. Summary of clinical evaluation and post-market clinical follow-up	10
5.1 Equivalence	10
5.2 Clinical investigations prior to CE mark	10
5.2.1 "Proof of concept"-study	11
5.2.1.1 Clinical investigation design	11
5.2.1.2 Results	13
5.2.1.3 Conclusion	14

5.2.2	Safety assessment study	15
5.2.2.1	Clinical investigation design	15
5.2.2.2	Results	17
5.2.2.3	Conclusion	18
5.2.3	Histopathological report	19
5.2.3.1	Clinical investigation design	19
5.2.3.2	Results	19
5.2.3.3	Conclusion	22
5.2.4	Prospective, open-label, clinical investigation of intra-articular polyacrylamide hydrogel injection in subjects with knee OA (“IDA”) (CIV ID: CIV-19-07-029146)	22
5.2.4.1	Clinical investigation design	23
5.2.4.2	Results	25
5.2.4.3	Conclusion	34
5.2.5	A multi-centre randomised, controlled, double-blind clinical investigation of intra-articular polyacrylamide hydrogel injection in subjects with knee osteoarthritis (“ROSA”) (CIV ID: CIV-19-01-026964)	34
5.2.5.1	Clinical investigation design	35
5.2.5.2	Results	37
5.2.5.3	Conclusion	46
5.3	Other clinical data	47
5.4	Summary of the clinical performance and safety	47
5.5	Post-market clinical follow-up	48
6.	Possible therapeutic alternatives	49
6.1	Intra-articular treatments	49
6.1.1	Corticosteroids	49
6.1.2	Hyaluronic acid	49
6.1.3	Noltrex	50
6.2	Surgery	50
7.	Suggested profile and training for users	51
8.	Harmonised standards and applied Common Specifications	51
9.	Revision history	54
10.	Bibliography	55
	Summary of Safety and Clinical Performance - Part II – Patients	61
1.	Device identification and general information	61
1.1	Device trade name	61
1.2	Manufacturer	61

1.3	Basic Unique Device Identification – Device Identifier (UDI-DI).....	61
1.4	Year when the device was first CE-marked	61
2.	Intended use of the device	61
2.1	Intended purpose, indication and intended patient group	61
2.2	Contraindications.....	62
3.	Device description	62
3.1	Device description and materials in contact with patient tissues	62
3.2	Information about medicinal substance, if any	62
3.3	Description of how the device is achieving its intended mode of action.....	62
3.4	Other devices or products intended to be used in combination with the device	63
4.	Risks and warnings.....	63
4.1	How potential risks have been controlled or managed	63
4.2	Remaining risks and undesirable effects	63
4.3	Warnings or precautions	64
5.	Summary of clinical evaluation and post-market clinical follow-up	64
5.1	Clinical background of the device	64
5.2	The clinical evidence for CE-marking	65
5.3	Safety	66
6.	Possible therapeutic alternatives	66
7.	Suggested profile and training for users	67

Abbreviations

Abbreviation	Explanation
ADE	Adverse device effect
AE	Adverse event
BOCF	Baseline observation carried forward
EMDN	European medical device nomenclature
ESCEO	European Society for Clinical and Economic Aspects of Osteoporosis, Osteoarthritis and Musculoskeletal Diseases
HA	Hyaluronic acid
IFU	Instructions for use
ITT	Intention to treat
MDR	Medical Device Regulation (EU) 2017/745
OA	Osteoarthritis
OMERACT-OARSI	Outcome Measures in Rheumatology - Osteoarthritis Research Society International
PAAG	Polyacrylamide hydrogel
PGA	Patient global assessment
PMCF	Post-market clinical follow-up
PMS	Post-market surveillance
QoL	Quality of life
ROM	Range of motion
SAE	Serious adverse event
SSCP	Summary of safety and clinical performance
TKA	Total knee arthroplasty
UKR	Unicompartmental knee replacement
WOMAC	Western Ontario and McMaster Universities Osteoarthritis Index

Summary of safety and clinical performance - Part I - Professional users

This summary of safety and clinical performance (SSCP) is intended to provide public access to an updated summary of the main aspects of the safety and clinical performance of Arthrosamid®.

The SSCP is not intended to replace the instructions for use (IFU) as the main document to ensure safe use of Arthrosamid®, nor is it intended to provide diagnostic or therapeutic suggestions to intended users or patients.

The following information is intended for users/healthcare professionals.

Following this information there is a summary intended for patients (Part II).

1. Device identification and general information

1.1 Device trade name

Arthrosamid®, catalogue no. 50020 (single pack) and 50021 (six-pack).

1.2 Manufacturer

Contura International A/S
Sydmarken 23
DK-2860 Soeborg
Denmark

1.3 Manufacturer single registration number (SRN)

DK-MF-000002018

1.4 Basic UDI-DI

5704101Arthrosamid3K

The Basic UDI-DI is the primary identifier of a device model.

1.5 Medical device nomenclature description

The European medical device nomenclature (EMDN) code for Arthrosamid® is P900403 “Non-resorbable filling and reconstruction devices”.

1.6 Class of device

Arthrosamid® is a class IIb implantable device according to Rule 8 set out in Annex VIII to the Medical Device Regulation (EU) 2017/745 (MDR).

1.7 Year when the first CE certificate was issued covering the device

Arthrosamid® was CE-marked under the Medical Device Directive 93/42/EEC in 2021.

1.8 Notified body name and identification number

The notified body conducting conformity assessment of Arthrosamid® is TÜV SÜD Product Service GmbH, with single identification number 0123.

2. Intended purpose, indication, target population and contraindications

2.1 Intended purpose, indication and target population

Arthrosamid® is intended to be used for symptomatic treatment of patients with knee osteoarthritis.

Arthrosamid® is intended for patients above the age of 18 years that have failed conservative treatments.

2.2 Contraindications and/or limitations

Arthrosamid® should not be injected:

- If an active skin disease or infection is present at or near the injection site.
- If the joint is infected or severely inflamed.
- If the patient has previously received treatment with a different non-absorbable injectable/implant in the same knee.
- If the patient has previously received treatment with an absorbable intra-articular injectable (such as hyaluronic acid) that is not yet expected to be absorbed, according to manufacturer's information for the specific product.
- If the patient has received a knee alloplasty or has any foreign material in the knee.
- If the patient has undergone knee arthroscopy within the last 6 months.
- In patients with uncontrolled haemophilia or in patients with uncontrolled anti-coagulation treatment.

Do not subsequently inject a different non-absorbable implant in the injected knee.

3. Device description

3.1 Description of device

Arthrosamid® is a non-absorbable, biocompatible, injectable, transparent hydrogel for intra-articular injection. Arthrosamid® consists of 2.5% cross-linked polyacrylamide and 97.5% non-pyrogenic water. Arthrosamid® is supplied in a pre-filled, sterile, 1 mL syringe, sealed with a Tip Cap in a Luer lock fitting.

Arthrosamid® is sterilised by moist heat. Intended for single use only.

Arthrosamid® is a permanent, non-biodegradable implant with an expected performance lifetime of 5 years.

3.1.1 Mode of action

- Arthrosamid® integrates into the synovial tissue of the inner joint capsule and decreases joint stiffness.
- Arthrosamid® forms a cushion-like membrane comprising the hydrogel incorporated in a matrix of fine tissue network covered by a synovial lining.
- Arthrosamid® diminishes pain and improves function of the knee affected by osteoarthritis.

The non-absorbable, non-bio-degradable and non-migratory characteristics of Arthrosamid® provide cushioning of the inner joint capsular tissue.

3.2 Previous generations of the device

There are no previous generations or variants of the device.

3.3 Any other devices or products intended to be used in combination with the device

Arthrosamid® is intended to be injected intra-articularly with a sterile 21G x 2 inches (0.8 x 50mm) needle with Luer lock fitting.

4. Residual risks and warnings

4.1 Residual risks and undesirable effects

Contura conducts risk management as an iterative process throughout the life cycle of Arthrosamid® to continuously manage all risks that arise in association with the device. A risk analysis has been performed for Arthrosamid®, in which 2 residual (after risk control) risks deriving from 1 reasonably foreseeable sequence or combination of events have been identified and evaluated as acceptable but significant. The harms associated with these 2 risks are transient symptoms not requiring professional

medical or surgical intervention which patients may experience during the first few weeks after the procedure (recovery period).

The harms associated with the residual risks are:

- Swelling/edema.
- Worsening of arthralgia.

Quantitative data on undesirable effects of Arthrosamid® are available from two clinical investigations, “IDA” (n=49) and “ROSA” (n=121) as well as from systematically gathered post-market surveillance (PMS) data up to and including 2023.DEC.31. A total of 5,098 patient packs (six-packs) were placed on the market per 2023.DEC.31.

The “IDA” and “ROSA” studies are described in sections 5.2.4 and 5.2.5, respectively.

Table 1 Frequencies of residual risks and undesirable effects in clinical trials and post-market

	“IDA” (n=49)	“ROSA” (n=121)	Post-market data (n=5,098)
Residual risks			
Swelling/edema	8%	11%	0.04%
Worsening of arthralgia	12%	16%	< 1/5,098*
Undesirable effects			
Injection-related pain	-	1%	0.08%
Swelling of Baker’s cyst	6%	-	< 1/5,098*
Infection	-	-	< 1/5,098*

*No product-related complications have been reported from post-market use.

4.2 Warnings and precautions

4.2.1 Warnings

- Do not inject Arthrosamid® intravascularly, extra-articularly or in the tissue of the capsule.
- Do not inject corticosteroids along with Arthrosamid®.
- Do not inject any pharmaceuticals or biological substances along with Arthrosamid®.

4.2.2 Precautions

Usual precautions associated with invasive joint procedures should be followed.

- Do not re-use Arthrosamid®. Re-use increases the risk of contamination and thereby increases the risk of infection.
- Do not re-sterilise Arthrosamid®.
- Do not use Arthrosamid® once expired.

- Do not inject Arthrosamid® into other parts of the body.
- Patients should avoid strenuous activity (e.g., tennis, jogging or long walks) during the first weeks after the injection.

Arthrosamid® should be used with caution in patients:

- with acute or chronic infections in other parts of the body.
- with autoimmune disorders.
- in immunosuppressive therapy.
- with a compromised immune system.
- suffering from uncontrolled diabetes.
- undergoing major dental work or an invasive medical procedure.

Safety and effectiveness have not been established for:

- patients under 18 years of age.
- pregnant or lactating women.
- patients with BMI ≥ 35 kg/m².
- single doses higher than 6 mL.

As with any invasive joint procedure, there is a risk of infection when injecting Arthrosamid® and therefore prophylactic antibiotics must be administered prior to injection.

4.3 Other information

Arthrosamid® has not been subject to any field safety corrective action.

5. Summary of clinical evaluation and post-market clinical follow-up

5.1 Equivalence

Conformity assessment of Arthrosamid® has not been based on equivalence to any other products.

5.2 Clinical investigations prior to CE mark

Clinical investigations of Arthrosamid® have been performed prior to CE-marking in accordance with the MDR.

A summary of the investigations is provided below.

5.2.1 “Proof of concept”-study

5.2.1.1 Clinical investigation design

This investigation was an investigator-initiated single-centre, open-label cohort “proof of concept” study. It was conducted in accordance with the declaration of Helsinki (WMA Declaration of Helsinki - Ethical Principles for Medical Research Involving Human Subjects) and published in 2018 [1].

Investigational device

Device: The “Hydrogel B” product Aquamid® Reconstruction was utilized for the study. Aquamid® Reconstruction is identical to Arthrosamid® but with another indication for use (Aquamid® Reconstruction was intended for soft tissue augmentation for reconstructive purposes, including corrections of facial lipoatrophy). In the following, the device is referred to as “Hydrogel B”.

Administration: The injection of “Hydrogel B” was performed under aseptic conditions and with local anaesthesia (inside the joint and directly at the injection site). A volume of 3 mL of “Hydrogel B” was injected under ultrasound guidance into the joint space using a 21G x 2” (0.8 x 50 mm) needle (3 x 1 mL was injected using the same needle). In subjects with insufficient effect, the injection was repeated after a minimum of 2 weeks, resulting in a total of 6 mL injected “Hydrogel B”.

Prophylactic antibiotics: A single oral dose of antibiotics (azithromycin 500 mg + moxifloxacin 400 mg) was administered 1-6 hours prior to injection.

Objectives and endpoints

The general study purpose was to obtain information on the effectiveness of “Hydrogel B” in subjects with osteoarthritis (OA) of the knee. The endpoints of the study were defined as:

Primary endpoint: To evaluate the change from baseline in WOMAC pain score at 4 months after injection of “Hydrogel B” into the knee joint. The WOMAC questionnaire is described below.

Secondary endpoints: To evaluate changes from baseline in the physical function and stiffness subscales of WOMAC as well as the total WOMAC score at 4, 7 and 13 months.

Effectiveness parameters

WOMAC score: The WOMAC questionnaire is a self-administered questionnaire consisting of 24 items, each scored on a 0-4 Likert scale. The 24 items are divided into three subscales: pain (5 items), stiffness (2 items) and physical function (17 items). A total WOMAC score is also calculated as the sum of all 24 items. The WOMAC scores (subscales and total) are converted to a 0–100 scale, with a higher score indicating a worse outcome.

Inclusion and exclusion criteria

The following inclusion and exclusion criteria were defined:

Inclusion criteria:

Adults with a clinical diagnosis of osteoarthritis of the knee according to the American College of Rheumatology.

Exclusion criteria:

- Active skin disease or infection at or near the injection site.
- Infected or severely inflamed joint.
- Previous treatment with degradable intra-articular injectable, such as hyaluronic acid, unless degradation of such product is complete.
- Previous treatment with a different non-absorbable injectable/implant.
- Subject has received a knee arthroplasty.
- Haemophilia patients.
- Subjects under the age of 18 years.

Study overview

The study comprises the following treatment and assessment time points:

- An initial treatment visit (Treatment visit 1).
- A second treatment visit approximately one month after Treatment visit 1 (Treatment visit 2).
- Follow-up visits at four months (Follow-up visit 1), seven months (Follow-up visit 2) and thirteen months (Follow-up visit 3). The scheduled visit windows were +/- 2 weeks for all follow-up visits.

A brief overview of the study procedure is outlined in [Table 2](#).

Table 2. “Proof of concept” study outline.

	Treatment visit 1	Treatment visit 2	Follow-up visit 1	Follow-up visit 2	Follow-up visit 3
Time point (months)	0	1	4	7	13
Injection of “Hydrogel B”	X	X*	-	-	-
WOMAC	X	-	X	X	X

*“Hydrogel B” injection performed if the effect of the initial dose was regarded as insufficient.

Study population

The study population consisted of subjects who were treated for knee OA and had valid records of pre-treatment primary effectiveness measure (WOMAC questionnaire). The subjects received the initial “Hydrogel B” injection at Treatment visit 1 and did not receive intra-articular corticosteroid

injection at any time point up to the 13 months study visit and did not receive intra-articular injection of “Hydrogel B” at any of the follow-up study visits (i.e., the 4- and 7-months follow-up visits).

The degree of OA was determined by the Kellgren-Lawrence (KL) classification, which is a grading from 0 to 4, with Grade 0 signifying no presence of OA and Grade 4 signifying severe OA [1]. The subjects in the study had the following KL grades: KL=1: 4%; KL=2: 24%; KL=3: 46%; and KL=4: 26%.

5.2.1.2 Results

Patient characteristics

Altogether, 122 knee OA patients received intra-articular treatment of 3 mL “Hydrogel B”, of which 38 did not complete a baseline WOMAC questionnaire, resulting in 84 subjects with pre-treatment WOMAC scores. The mean age at the first “Hydrogel B” injection was 68 years (range 36-91). 48 (57%) of the patients were female. Baseline WOMAC scores are given in [Table 3](#).

Table 3. WOMAC scores at baseline.

	Mean (SD)	Range
WOMAC pain	44.3 (17.2)	5.0 – 75.0
WOMAC stiffness	44.9 (22.0)	0.0 – 87.5
WOMAC function	42.2 (20.7)	0.0 – 88.2
WOMAC total	42.8 (18.6)	0.0 – 84.4

Treatment and follow-up

Of the 84 subjects with pre-treatment WOMAC scores, 75 subjects attended the second treatment visit, of which 72 subjects received a second intra-articular injection of 3 mL of “Hydrogel B” (3 subjects did not attend the second treatment visit, but continued participation). A cumulative effect of “Hydrogel B” has therefore been studied indirectly.

WOMAC follow-up data were available for 62 subjects at 4 months, for 59 subjects at 7 months and for 56 subjects at 13 months.

Effectiveness parameters

WOMAC scores

[Table 4](#) summarises the WOMAC sub-scores (pain, stiffness, function) and the total WOMAC score. At the 4 months follow-up, there was a statistically significant decrease in the WOMAC pain sub-score of 14.6 points (95% CI -18.9 to -10.2) ($p < 0.001$). The improvement in pain sub-score persisted over time with a mean reduction of 16.0 and 15.7 at 7 and 13 months, respectively. Furthermore, WOMAC stiffness sub-score, WOMAC function sub-score and the total WOMAC score were all statistically significantly reduced at all time points (4, 7 and 13 months).

Table 4. Changes from baseline in WOMAC pain, stiffness and function sub-scores, and in WOMAC total score (n=84).

	4 months (n=62)		7 months (n=59)		13 months (n=56)	
	Mean (95% CI)	P-value	Mean (95% CI)	P-value	Mean (95% CI)	P-value
WOMAC pain	-14.6 (-18.9 to -10.2)	<.0001	-16.0 (-20.4 to -11.6)	<.0001	-15.7 (-20.2 to -11.2)	<.0001
WOMAC stiffness	-12.3 (-17.7 to -6.9)	<.0001	-13.3 (-18.8 to -7.8)	<.0001	-12.1 (-17.8 to -6.4)	<.0001
WOMAC function	-13.1 (-17.4 to -8.7)	<.0001	-12.2 (-16.6 to -7.8)	<.0001	-9.4 (-14.0 to -4.9)	<.0001
WOMAC total	-13.4 (-17.5 to -9.2)	<.0001	-13.1 (-17.3 to -8.8)	<.0001	-10.9 (-15.2 to -6.6)	<.0001

Achievement of study endpoints

The study results as related to the study endpoints are summarised in [Table 5](#).

Table 5. Study results as related to study endpoints.

Study endpoints		Study results
Primary endpoint	To evaluate the change from baseline in pain (WOMAC) at 4 months after injection of “Hydrogel B” into the knee joint.	A statistically significant improvement in pain at 4 months after injection of “Hydrogel B” into the knee joint was demonstrated ($p<0.0001$).
Secondary endpoints	To evaluate changes from baseline in the physical function and stiffness sub-scales of WOMAC as well as the total WOMAC score at 4, 7, and 13 months.	Improvement was seen in the other WOMAC domains (stiffness, function, and total) after 4 months. The effect seems to be long-term with effectiveness maintained for 13 months.

5.2.1.3 Conclusion

This clinical study on Arthrosamid® was a single-centre, open-label cohort “proof of concept” study including 84 subjects with pre-treatment WOMAC scores. Evaluation of the primary endpoint revealed that the WOMAC pain sub-score was significantly improved four months after the initial “Hydrogel B” injection as compared to baseline. Similarly, the WOMAC total score and the stiffness and physical function sub-scores (secondary endpoints) were all improved. The effect appeared to persist over time, up to the last observation at 13 months.

5.2.2 Safety assessment study

5.2.2.1 Clinical investigation design

General information

This was an investigator-initiated retrospective cohort study set up to assess the safety of Arthrosamid® (“Hydrogel B”) in patients with knee OA. The study was published in 2019 [2].

Investigational device

Device: The “Hydrogel B” product Aquamid® Reconstruction was utilized for the study. Aquamid® Reconstruction is identical to Arthrosamid® but with another indication for use (Aquamid® Reconstruction was designed for soft tissue augmentation for reconstructive purposes, including corrections of facial lipoatrophy). In the following, the device is referred to as “Hydrogel B”.

Administration: The injection of “Hydrogel B” was performed under aseptic conditions and with local anaesthesia (inside the joint and directly at the injection site). A volume of 3 mL of “Hydrogel B” was injected under ultrasound guidance into the joint space using a 21G x 2” (0.8 x 50 mm) needle (3 x 1 mL was injected using the same needle). In subjects with insufficient effect, the injection was repeated after a minimum of 2 weeks, resulting in a total of 6 mL injected “Hydrogel B”.

Prophylactic antibiotics: A single oral dose of antibiotics (azithromycin 500 mg + moxifloxacin 400 mg) was administered 1-6 hours prior to injection.

Study objective

The overall study purpose was to obtain information on the safety of “Hydrogel B” in subjects with OA of the knee.

Safety parameters

Adverse events: An adverse event is defined as any untoward medical occurrence or any untoward clinical signs (including an abnormal laboratory finding) in patients whether or not related to the investigational medical device. Adverse events were graded as mild (transient and easily tolerated), moderate (causing discomfort and interrupting normal activities) or severe (incapacitating, life-threatening or causing considerable interference with normal activities).

Serious adverse events: A serious adverse event is an adverse event that:

- led to death, injury or permanent impairment to a body structure or body function.
- led to a serious deterioration in health of the subject, that either resulted in:
 - a life-threatening illness or injury, or
 - a permanent impairment of a body structure or body function, or
 - in-subject hospitalisation or prolongation of existing hospitalisation, or
 - in medical or surgical intervention to prevent life-threatening illness.
- led to foetal distress, foetal death or congenital abnormality or birth defect.

Effectiveness parameter

Treatment effect questionnaire: To assess the subject's perceived effect of the treatment, the subjects were given three options for describing the current status of their knee compared to the pre-treatment status: unchanged / worse / better. If the answer was "worse" the subject was asked to perform a rating according to a 7-point scale where -1 denotes "almost the same, hardly any worse at all" and -7 denotes "a very great deal worse". If the answer was "better" the subject was asked to perform a rating according to a 7-point scale where 1 denotes "almost the same, hardly any better at all" and 7 denotes "a very great deal better".

The subject's perceived treatment effect was evaluated separately for pain, physical function and global perceived effect.

Study overview

All subjects had been treated with "Hydrogel B" for knee OA between 2010 and 2017. After informed consent was given, subjects underwent a medical examination and medical records were gathered. Subjects who declined the clinical examination were invited to be interviewed by telephone, and detailed descriptions of any adverse events were recorded. Thus, the study comprised one single hospital visit / telephone interview per subject.

Subject interview: All subjects were asked to describe any pain and mobility limitations, skin conditions, systemic events (not knee-related events) that they experienced in relation to or since the treatment with "Hydrogel B". If any adverse events were described, the subject was asked for duration and severity.

Medical records review / surgical reports: The investigator assessed any adverse event in detail including the approximate date of onset, description, severity, duration and outcome, relationship of the adverse event to investigational device, and any action(s) taken.

For the subset of subjects that had received knee replacement surgery since receiving treatment with "Hydrogel B", the surgery reports were reviewed for descriptions of abnormal macroscopic findings. In particular it was investigated if the presence of "Hydrogel B" was noticed and if there were signs of fibrosis or unidentified foreign bodies.

Clinical examination: A clinical examination of the participants' knees was performed to assess if abnormalities not otherwise found in relation to knee OA were present.

The clinical examination was specifically focused on the potential development of intra-articular fibrosis. Thus, an evaluation of knee joint range of motion was performed. In the case of a clinically significant reduction in flexion range of motion, it was judged and recorded whether this was expected in the natural course of the disease (knee OA) and if it was caused by either bony or soft tissue changes (fibrosis).

5.2.2.2 Results

A total of 128 subjects were referred to the OA outpatient clinic for potential inclusion in the study. Of these, 91 responded to the invitation and represent the retrospective safety population. A majority of the subjects visited the clinic (n=74), whereas 17 subjects participated by way of telephone interview.

The mean time since the initial “Hydrogel B” injection was 24 months (range: 4 to 87). [Table 6](#) shows the patient characteristics at the time of the study.

Table 6. Patient characteristics.

Number of subjects (n)	91
Age (years)	Mean: 70 Range: 38 - 94
Gender (n)	Female: 46 (51%) Male: 45 (49%)
Time since first “Hydrogel B” injection (months)	Mean: 24 Range: 4 - 87

Adverse events

In this retrospective assessment, 25 subjects (27%) reported that they had experienced one or more adverse events or discomfort in relation to or since their treatment with “Hydrogel B”. The total number of adverse events reported by these 25 subjects was 41. An overview of the adverse events as to frequency, severity and duration is given in [Table 7](#). The most frequent adverse events were a sensation of distension (15 subjects) or pain from the target knee (7 subjects).

Most of the subjects who had a sensation of distension described it as mild or moderate. Three subjects rated the sensation of distention as severe. All but one of the subjects described the sensation as short-lived (days to weeks). Pain in the knee was described as severe by 4 subjects. In none of the subjects the pain persisted for a prolonged period.

The majority of the patient-reported adverse events were judged by the investigator to be normal symptoms of knee OA.

Two subjects reported to have sought medical assistance at the treating clinic because of knee pain after the initial treatment. One was treated with analgesics and was in remission after a week. The other had intra-articular effusion and was treated with recurrent arthrocentesis, one and three weeks after the initial treatment, respectively. No lasting reactions and no allergic reactions or intra-articular infections were reported.

Table 7. Patient-reported adverse events.

Adverse event	n (%)	Severity (n)			Duration (n)		
		Mild	Moderate	Severe	Days	Weeks	Months
Pain from target knee	7 (17.1%)	2	1	4	3	4	0
Other pain from extremity	1 (2.4%)	0	0	1	1	0	0
Muscle discomfort/pain	2 (4.9%)	2	0	0	2	0	0
Soreness	3 (7.3%)	3	0	0	2	1	0
Burning sensation	1 (2.4%)	1	0	0	1	0	0
Sensation of distension	15 (36.6%)	10	2	3	6	8	1
Skin or joint pricking sensation	3 (7.3%)	1	1	1	1	1	1
Numbness	1 (2.4%)	1	0	0	1	0	0
Cold sensation	1 (2.4%)	1	0	0	0	0	1
Heat sensation	1 (2.4%)	1	0	0	0	1	0
Reduced range of motion	4 (9.8%)	1	2	1	1	2	1
Stiffness	2 (4.9%)	2	0	0	0	1	1

Serious adverse events

No serious adverse events were reported.

Clinical examination

From the clinical examinations of the 74 subjects who visited the osteoarthritis outpatient clinic, the investigator found no signs of reduced knee range of motion (ROM) other than what is expected in a population of knee OA patients.

The investigator found no indication of fibrosis formation significant enough to be detected through a clinical ROM examination in any of the assessed subjects.

Surgical reports

Fifteen of the subjects had undergone arthroplasty surgery for knee OA of the treated knee subsequent to the “Hydrogel B” injections. The average time from treatment with “Hydrogel B” to surgery was 19 months (range: 5 to 58 months). From the pre-operative examinations there were no indications of knee fibrosis in any of the subjects. There were no descriptions of abnormalities that could be associated with prior treatment with “Hydrogel B” in the surgical reports. None of the reports mentioned the presence of “Hydrogel B”, unidentified foreign bodies or macroscopical signs of fibrosis.

5.2.2.3 Conclusion

This clinical study on Arthrosamid® (“Hydrogel B”) was a retrospective cohort study including 91 subjects. Evaluation of safety data showed no serious adverse events related to “Hydrogel B” injection. No allergic reactions or intra-articular infections were observed.

Most of the adverse events reported by the subjects can be described as normal symptoms of knee OA and were thus most probably not related to “Hydrogel B”. An exception is the sensation of distention after treatment which was felt by 15 subjects. Most (14 subjects, corresponding to 93%) describe this sensation as short-lived (days to weeks). The cause of this adverse event is not clear.

There were no indications of fibrosis after “Hydrogel B” injection.

5.2.3 Histopathological report

This investigation is described in detail in the report: “Case series: Report on a prospective histopathological study on biopsies from seven total arthroplasty knees (TKA) treated with polyacrylamide hydrogel (PAAG)” [3].

5.2.3.1 Clinical investigation design

An investigator-initiated prospective histopathological study was performed on tissue removed from subjects during total knee arthroplasty (TKA). The subjects had previously received treatment with “Hydrogel B”. The surgical procedure at TKA includes resection of diseased or damaged bone and cartilage from the knee, followed by attachment of metal and polyethylene prostheses.

Upon resection, soft tissue, cartilage and decalcified bone were cut into smaller sections and processed for histological examination. Haematoxylin-eosin-stained sections were examined.

The augmentation effect of “Hydrogel B” was estimated as a combination of increased average thickness of the synovial membrane and extent of “Hydrogel B” integrated foci and scored on a scale between 1 and 3.

The integration grade was judged as a combination of network mask-size and number of cells within this network and scored on a scale between 1 and 3.

5.2.3.2 Results

A total of seven subjects were included in the study (Table 8). Subject age ranged between 35 and 84 years, and all but one of the subjects were female. In this study, the stages of OA are given as OA grade¹ (determined by the surgeon) and/or KL score. Both are graded from 0 to 4, with Grade 0 signifying no presence of OA and Grade 4 signifying severe OA. OA was the indication for “Hydrogel B” injection, and subsequent total knee arthroplasty, in six of the subjects. The seventh subject had suffered from a knee injury at the age of 13. “Hydrogel B” had been injected twice in a majority of the subjects over a mean period of 16 months (range 5-33) prior to knee arthroplasty. The reason for arthroplasty was, in most cases, recurrence of symptoms.

The histological findings are summarised in Table 9. The augmentation effect of “Hydrogel B” was estimated as a combination of increased average thickness of the synovial membrane and the extent of “Hydrogel B” integrated foci, and scored on a scale between 1 and 3, with 3 being the highest augmentation effect. The integration grade was judged as a combination of network mask-size and

¹ <https://www.healthline.com/health/osteoarthritis-stages-of-oa-of-the-knee>

number of cells within this network, and scored on a scale between 1 and 3, with 3 being the highest integration grade.

Table 8. Patient characteristics.

Patient	Age / Sex	Disease	OA grade	KL score*	Previous treatment	"Hydrogel B" volume	Time since "Hydrogel B" injection (months)
Case 1	35 years - female	Injury	1	-	Surgery	5 + 3 mL	18
Case 2	84 years - female	OA	2	3	Steroid	3 + 3 mL	5
Case 3	77 years - female	OA	2	-	None	5 + 5 mL	29
Case 4	69 years - female	OA	2	2	None	3 mL	33
Case 5	67 years - female	OA	3	3	None	3 + 3 mL	7
Case 6	74 years - female	OA	3	2	None	3 + 3 mL	5
Case 7	82 years - male	OA	2	3	Steroid	3 + 3 mL	16

*KL score was not available in case 1 and case 3.

Table 9. Summary of histological findings.

Subject	Biopsies with gel	Augmentation effect of gel (1-3)	Integration grade (1-3)	Summary of histological findings
Case 1	3/16	1	2	Scattered small gel deposits just beneath the synovial lining. Signs of OA with cartilage fragments, synovial hyperplasia, diffuse increase in inflammatory cells of the synovial membrane.
Case 2	6/25	1-2	1	A few gel deposits just beneath the synovial lining. A few spots of hyalinised connective tissue, a few spots of inflammatory cells in the synovial membrane. No areas of sclerosing / scar-forming fibrosis or necrosis.
Case 3	38/48	3	3	Large quantities of fully integrated gel within the synovial membrane, within the connective tissue and beneath the synovial membrane. Foci of inflammatory cells and streaks of connective tissue fibres.
Case 4	12/26	3	3	Fully integrated gel with scattered inflammatory cells and fibroblasts.
Case 5	3/15	1-2	3	Irregularly folded capsule tissue in areas with large papillary folds. Integrated gel was seen sub-synovially.
Case 6	10/11	1	2	Well integrated gel within the synovial membrane and beneath the synovial lining. Diffuse infiltration of inflammatory cells in the connective tissue.
Case 7	11/18	2	2	Large condensations of inflammatory cells around the gel deposits. Small foci with bleeding and necrosis.

In essence, a similar histological pattern was seen in all seven cases. "Hydrogel B" was found to be integrated in the synovial membrane, and outer synovial lining cells had entered the gel and established a de novo lining layer. The integrated gel appears to contribute with a considerable augmentation effect to the inner side of the capsule, facing the cavity. Thus, extensive papillary folding of the synovial membrane was seen in all cases except one (case 1). In one knee (case 1) the

gel was found immediately adjacent to both vital and necrotic bone trabecula, where it might act as an extra padding.

None of the seven cases showed evidence of firm fibrosis or scar formation. The integration of the gel resulted in a network of very fine fibres of tissue network, but none of the many biopsies showed any signs of shrinkage or contraction within or surrounding the gel deposits that might give indication of fibrosis or scarring, which are well-known parameters in OA.

Examples of the histological findings are shown in [Figure 1](#) and [Figure 2](#).


Figure 1. Case 3. Overview of part of the synovial membrane showing papillary proliferations augmented with “Hydrogel B” (I). Note the hyperplastic synovial lining (circled). Slender bone trabecula (B/bone) and intervening marrow spaces augmented with integrated gel (G) (II and III). Necrotic and regenerating cartilage (RG) or cartilage/bone (C) with integrated gel (G) lying in adjacent marrow spaces (IV).


Figure 2. Case 7. Synovial papillary projections with integrated gel of different volumes and frequencies within the synovial membrane (arrows) (I and II). The three black circles illustrate different contents of gel within the different synovial projections. The smallest circle displays no gel, the middle size one contains some gel, and the large circle contains gel in all tissue projections (I and II). Condensations of chronic inflammatory cells are seen within the red circle (I). Folded rim of integrated gel is shown in (III) and a close-up photo of an area with gel integration and chronic inflammatory cells in surrounding connective tissue is seen in (IV).

5.2.3.3 Conclusion

This pathohistological study included knee biopsies from seven subjects, previously treated with “Hydrogel B”, obtained in conjunction with total knee arthroplasty. The overall finding of the investigation is that “Hydrogel B” was found to be integrated in the synovial membrane and covered by synovial surface cells, irrespective of the time point after “Hydrogel B” injection (5 – 33 months, mean: 16 months). In addition, extensive papillary folding of the synovial membrane was observed. No evidence of firm fibrosis or scar formation was seen.

5.2.4 Prospective, open-label, clinical investigation of intra-articular polyacrylamide hydrogel injection in subjects with knee OA (“IDA”) (CIV ID: CIV-19-07-029146)

Contura is sponsoring a multi-centre, prospective, open-label clinical study evaluating the effectiveness and safety of Arthrosamid® in subjects with knee OA. The investigation consists of a main study followed by an extension study. 49 subjects were enrolled over a period of one month (August-September 2019). The study was conducted according to ISO 14155 and all applicable regulatory requirements, at three clinical sites in Denmark, and prior to initiation of the investigation, approvals were obtained from the Danish Medicines Agency and the Independent Ethics Committee.

Subjects have been followed for 4 years. Data evaluated at the 6- and 12-month time point were published in 2021 [4,5], 2-year data were published in 2022 [6], and 3-year data were published in 2023 [7,8].

5.2.4.1 Clinical investigation design

Clinical investigation objectives and endpoints

The primary objective of the main study was to evaluate the effectiveness of one injection of 6 mL (6 x 1 mL) Arthrosamid® on pain over 3 months (13 weeks) in subjects with mild to severe knee OA (KL grade 2-4).

The primary endpoint for the main study was:

- Change in the 0-100 points normalised WOMAC pain subscale from baseline to week 13.

The secondary objectives were to evaluate the effectiveness of 6 mL Arthrosamid® on knee OA-related subject reported symptoms, function during activities of daily living, and patient global assessment (PGA) of impact of OA over 1 month (4 weeks) and 3 months (13 weeks).

The secondary endpoints for the main study were:

Changes from baseline over 1 month (4 weeks) and 3 months (13 weeks) in:

- WOMAC stiffness and physical function-related subscales.
- PGA of impact of OA.
- Outcome Measures in Rheumatology - Osteoarthritis Research Society International (OMERACT-OARSI) response.
- WOMAC pain-related subscale (1 month).

The extension study only includes secondary endpoints. These are changes from baseline to weeks 26 and 52, year 2, 3, 4, and 5 in:

- WOMAC pain, stiffness, and physical function-related subscales.
- PGA of impact of OA.
- OMERACT-OARSI response (only week 26 and 52).

The safety objectives (for both main and extension study) are to evaluate safety of Arthrosamid® in subjects with knee OA.

Safety endpoints (for both main and extension study) were:

- Incidence of adverse events (AEs).
- Incidence of adverse device effects (ADEs).

Study treatment

48 of the 49 subjects received the planned 6 mL injection of the device. Due to a device deficiency one subject only received approximately 5 mL of the device. This subject completed year 4 of the extension study.

All subjects were given a local anaesthetic prior to the injection, and all received prophylactic treatment with azithromycin and moxifloxacin.

Patient population

Main inclusion criteria were:

- Signed informed consent form.
- >18 years old.
- Clinical diagnosis of knee OA according to American College of Rheumatology criteria and a KL score of 2, 3 or 4.

Main exclusion criteria were according to the contraindications for Arthrosamid®:

- If an active skin disease or infection was present at or near the injection site.
- If the joint was infected or severely inflamed.
- If a degradable intra-articular injectable such as hyaluronic acid had been used within the last 6 months.
- If the subject had previously received treatment with a different non-absorbable injectable/implant.
- If the subject had received a knee arthroplasty.
- Haemophilia patients.

In order to adequately distinguish the response to Arthrosamid® from any background “noise” and to separate out any confounding factors, additional exclusion criteria (other than what is already taken into account with warnings and precautions) were implemented:

- Previous intra-articular injection of polyacrylamide gel in the target knee.
- Significant valgus/varus deformity of the knee, ligamentous laxity or meniscal instability.
- Diseases in target knee other than OA.
- Intra-articular injection of any substance other than hyaluronic acid (e.g., corticosteroids) in the target knee within the last 3 months.
- History of surgery in the target knee within the past 6 months.
- Symptomatic OA of the hips, spine or ankle, that interferes with the evaluation of the target knee.
- Planned surgery on any lower extremity.
- Clinically significant venous or lymphatic stasis present in the legs.
- Any unstable or severe cardio-vascular disease.
- Any other contraindication to intra-articular injection.

- Any significant medical condition (e.g., significant psychiatric or neurological disorders or active alcohol/drug abuse) that is unstable/poorly controlled or other factors that may interfere with study participation.
- Treatment with systemic steroids with daily doses equivalent of >7.5 mg prednisolone.
- Significant change in physiotherapy in lower extremities related to OA within the previous month.
- Fibromyalgia.
- Inflammatory or other disease/condition which may affect the knee joint (e.g., rheumatoid arthritis, psoriatic arthritis and ankylosing spondylitis among others).
- Any other condition that in the opinion of the investigator puts a potential participant at risk or otherwise precludes participation in the investigation.
- Participation in any experimental device study within 6 months prior to enrolment, or participation in an experimental drug study within 1 month prior to enrolment.

5.2.4.2 Results

Disposition of subjects

50 subjects were screened, and 49 subjects were enrolled into the study (intention to treat [ITT] analysis set). All subjects completed the main study and entered the 1-year extension phase. 46 subjects completed the 1-year extension.

A total of 35 subjects entered year 2 of the extension study, and 32 subjects completed the 2-year extension. One subject withdrew consent (withdrawal by subject) during year 2, and 2 subjects were withdrawn with reason for withdrawal reported as 'other'. Of the 3 subjects withdrawn from year 2 of the extension, 2 had target knee replacement on the date of withdrawal. For the remaining subject, no further information is available regarding the reason for withdrawal.

A total of 29 subjects completed the 3-year extension (visit 8). One subject withdrew consent (withdrawal by subject) during year 3, and 2 subjects were withdrawn due to knee replacement surgery.

A total of 27 subjects completed the 4-year extension (visit 9). Two subjects withdrew consent (withdrawal by subject) during year 4.

Effectiveness endpoints

WOMAC pain subscale

The scores for the transformed WOMAC pain subscale at baseline, week 13, week 26, week 52, year 2, year 3 and year 4 are summarised for the ITT analysis set together with the respective changes from baseline in [Table 10](#). The changes from baseline are shown graphically in [Figure 3](#) for the ITT analysis set, both as observed and using the "baseline observation carried forward" (BOCF) approach, in which any missing values at 2 years were replaced by the participants' baseline values, and for the year-2 extension subset.

The analysis of the change from baseline in the transformed WOMAC pain subscale (0-100) showed a significant reduction ($p < 0.0001$) of 22.0 points from baseline to year 4 for the full ITT analysis set

and 21.8 points for the year-2 extension subset, indicating reduced pain in the treated knee. The most conservative approach, BOCF, also resulted in a significant reduction ($p < 0.0001$) of 13.0 points from baseline to year 4.


Figure 3. Mean plot of change from baseline (CfB) in transformed WOMAC pain subscale (0-100). Error bars represent standard deviation.

Table 10. Summary of transformed WOMAC pain subscale (0-100).

	Actual values	Change from baseline
ITT analysis set (N,%)	49 (100.0)	
Baseline		
N	49	
Mean (SD)	50.3 (11.8)	
Median	50.0	
Min - Max	20 - 75	
Week 13		
N	48	48
Mean (SD)	32.1 (18.2)	-18.0 (18.5)
Median	35.0	-20.0
Min - Max	0 - 75	-55 - 35
Week 26		
N	46	46
Mean (SD)	28.9 (18.6)	-21.7 (19.6)
Median	30.0	-22.5
Min - Max	0 - 70	-60 - 25
Week 52		
N	45	45
Mean (SD)	32.1 (20.0)	-18.4 (17.7)
Median	30.0	-20.0
Min - Max	0 - 70	-60 - 20
Year 2		
N	32	32
Mean (SD)	28.9 (19.0)	-20.3 (18.8)
Median	25.0	-22.5
Min - Max	0 - 70	-50 - 15
Year 3		
N	29	29
Mean (SD)	28.8 (18.0)	-20.5 (17.7)
Median	30.0	-15.0
Min - Max	0 - 70	-60 - 5
Year 4		
N	27	27
Mean (SD)	24.6 (17.0)	-23.5 (21.4)
Median	25.0	-25.0
Min - Max	0 - 55	-65 - 20

N: Number of subjects, SD: Standard deviation.

WOMAC stiffness subscale

The scores for the transformed WOMAC stiffness subscale at baseline, week 13, week 26, week 52, year 2, year 3 and year 4 are summarised for the ITT analysis set together with the respective changes from baseline in [Table 11](#). The mean changes from baseline are shown graphically in [Figure](#)

4 for the ITT analysis set, both as observed and using the BOCF approach, and for the year-2 extension subset.

The analysis of the change from baseline in the transformed WOMAC stiffness subscale (0-100) showed a significant reduction ($p < 0.0001$) of 23.5 points from baseline to year 4 for the full ITT analysis set and 22.5 points for the year-2 extension subset, indicating reduced stiffness in the treated knee. The most conservative approach, BOCF, also resulted in a significant reduction ($p < 0.0001$) of 13.5 points from baseline to year 4.


Figure 4. Mean plot of change from baseline (CfB) in transformed WOMAC stiffness subscale (0-100). Error bars represent standard deviation.

Table 11. Summary of transformed WOMAC stiffness subscale (0-100).

	Actual values	Change from baseline
ITT analysis set (N,%)	49 (100.0)	
Baseline		
N	49	
Mean (SD)	55.6 (17.5)	
Median	62.5	
Min - Max	0 - 88	
Week 13		
N	48	48
Mean (SD)	34.6 (21.0)	-20.8 (19.9)
Median	25.0	-25.0
Min - Max	0 - 75	-63 - 38
Week 26		
N	46	46
Mean (SD)	37.2 (21.5)	-17.9 (21.2)
Median	31.3	-18.8
Min - Max	0 - 75	-50 - 50
Week 52		
N	45	45
Mean (SD)	43.9 (24.8)	-11.7 (20.9)
Median	50.0	-12.5
Min - Max	0 - 88	-50 - 50
Year 2		
N	32	32
Mean (SD)	35.9 (23.7)	-19.1 (19.8)
Median	37.5	-12.5
Min - Max	0 - 75	-75 - 13
Year 3		
N	29	29
Mean (SD)	37.1 (19.9)	-16.8 (15.4)
Median	37.5	-12.5
Min - Max	0 - 75	-50 - 13
Year 4		
N	27	27
Mean (SD)	28.7 (24.2)	-24.5 (22.1)
Median	25.0	-25.0
Min - Max	0 - 75	-63 - 25

N: Number of subjects, SD: Standard deviation.

WOMAC physical function subscale

The scores for the transformed WOMAC physical function subscale at baseline, week 13, week 26 and week 52, year 2, year 3 and year 4 are summarised for the ITT analysis set together with the respective changes from baseline in [Table 12](#). The changes from baseline are shown graphically in

Figure 5 for the ITT analysis set, both as observed and using the BOCF approach, and for the year-2 extension subset.

The analysis of the change from baseline in the transformed WOMAC physical function subscale (0-100) showed a significant reduction ($p < 0.0001$) of 17.5 points from baseline to year 4 for the full ITT analysis set and 17.2 points for the year-2 extension subset, indicating an improved physical function in the treated knee. The most conservative approach, BOCF, also resulted in a significant reduction ($p < 0.0001$) of 9.7 points from baseline to year 4.


Figure 5. Mean plot of change from baseline (CfB) in transformed WOMAC physical function subscale (0-100). Error bars represent standard deviation.

Table 12. Summary of transformed WOMAC physical function subscale (0-100).

	Actual values	Change from baseline
ITT analysis set (N,%)	49 (100.0)	
Baseline		
N	49	
Mean (SD)	46.6 (16.1)	
Median	45.6	
Min - Max	9 - 87	
Week 13		
N	48	48
Mean (SD)	29.4 (18.0)	-17.1 (16.1)
Median	28.7	-19.8
Min - Max	0 - 60	-44 - 25
Week 26		
N	46	46
Mean (SD)	28.6 (20.4)	-18.4 (17.0)
Median	28.7	-20.2
Min - Max	0 - 62	-47 - 21
Week 52		
N	43	43
Mean (SD)	31.0 (22.3)	-16.0 (13.5)
Median	26.5	-16.2
Min - Max	0 - 69	-40 - 16
Year 2		
N	32	32
Mean (SD)	23.9 (18.0)	-18.8 (16.0)
Median	22.1	-17.6
Min - Max	0 - 62	-50 - 13
Year 3		
N	29	29
Mean (SD)	26.9 (18.8)	-15.2 (18.1)
Median	25.0	-11.8
Min - Max	0 - 62	-54 - 13
Year 4		
N	27	27
Mean (SD)	23.9 (19.1)	-17.6 (21.6)
Median	19.1	-17.6
Min - Max	0 - 62	-69 - 26

N: Number of subjects, SD: Standard deviation.

Patient global assessment

The PGA was in response to the question “How much does the knee osteoarthritis (target knee) as a whole affect your life at present?”. The response was indicated on a 100 mm visual analogue scale between “Not at all” (0) and “The worst possible” (100). The scores for the PGA at baseline, week 13,

Week 26, week 52, year 2, year 3 and year 4 are summarised for the ITT analysis set together with the respective changes from baseline in [Table 13](#). The changes from baseline are shown graphically in [Figure 6](#) for the ITT analysis set, both as observed and using the BOCF approach, and for the year-2 extension subset.

The analysis of the change from baseline in the PGA showed a significant reduction of 21.2 points ($p < 0.0001$) from baseline to year 4 for the full ITT analysis set and 19.1 points ($p = 0.0001$) for the year-2 extension subset, indicating a reduction in the mean impact of knee OA on the subjects' lives. The most conservative approach, BOCF, also resulted in a significant reduction ($p < 0.0001$) of 11.9 points from baseline to year 4.


Figure 6. Mean patient global assessment. Error bars represent standard deviation.

Table 13. Summary of patient global assessment.

	Actual values	Change from baseline
ITT analysis set (N,%)	49 (100.0)	
Baseline		
N	49	
Mean (SD)	61.1 (18.3)	
Median	65.0	
Min - Max	22 - 100	
Week 13		
N	48	48
Mean (SD)	42.6 (24.7)	-18.6 (24.4)
Median	40.0	-19.5
Min - Max	2 - 95	-67 - 53
Week 26		
N	44	44
Mean (SD)	43.9 (27.6)	-17.2 (27.5)
Median	41.0	-16.0
Min - Max	0 - 97	-89 - 46
Week 52		
N	45	45
Mean (SD)	43.3 (27.1)	-17.7 (22.7)
Median	36.0	-15.0
Min - Max	1 - 95	-67 - 24
Year 2		
N	32	32
Mean (SD)	33.9 (23.6)	-21.4 (22.5)
Median	32.5	-20.5
Min - Max	0 - 84	-78 - 31
Year 3		
N	29	29
Mean (SD)	37.2 (28.0)	-18.3 (30.5)
Median	27.0	-18.0
Min - Max	0 - 89	-88 - 38
Year 4		
N	27	27
Mean (SD)	33.2 (19.4)	-21.6 (22.2)
Median	28.0	-20.0
Min - Max	7 - 68	-66 - 27

N: Number of subjects, SD: Standard deviation.

OMERACT-OARSI response

The OMERACT-OARSI criteria is a set of responder criteria, which present the results of changes after treatment of OA in three symptomatic domains (pain, function, and PGA) as a single variable. Change in OMERACT-OARSI response was assessed. More than half of the subjects were positive

responders to the OMERACT-OARSI criteria at week 4 (53.1%), week 13 (64.6%), week 26 (67.4%), and week 52 (62.2%). The OMERACT-OARSI response was only evaluated up to week 52.

Safety endpoints

Adverse events

Thirty-six subjects reported a total of 80 AEs up to and including the 4-year extension visit. Nine of the events were reported as serious adverse events (SAEs). No deaths were reported, and one of the AEs (arthralgia, unlikely related) led to withdrawal from the study.

The most frequently reported AEs were arthralgia (18 events in 15 subjects), corona virus infection (7 events in 7 subjects), joint swelling (6 events in 6 subjects), back pain (4 events in 4 subjects) and synovial cyst (4 events in 4 subjects). Almost all the AEs were mild (47 events in 30 subjects) or moderate (25 events in 18 subjects) in intensity. Eight severe events were reported in 7 subjects: arthralgia, cerebrovascular accident, corona virus infection, lumbar spinal stenosis, presyncope, prostate cancer uterine prolapse and complicated appendicitis.

Adverse device effects

An adverse device effect (ADE) is defined as an AE related to the use of an investigational medical device.

Of the 80 AEs, 14 ADEs occurred in 12 subjects. No serious ADEs and no unanticipated serious ADEs occurred. The most frequently reported ADEs were arthralgia (6 events in 6 subjects), joint swelling (4 events in 4 subjects) and synovial cyst (3 events in 3 subjects). No other event was reported in more than 1 subject. Two events of arthralgia in 2 subjects had moderate intensity, while all other ADEs were mild.

5.2.4.3 Conclusion

The investigation provides valid scientific evidence that Arthrosamid® remains safe and effective for its intended use 4 years after injection.

5.2.5 A multi-centre randomised, controlled, double-blind clinical investigation of intra-articular polyacrylamide hydrogel injection in subjects with knee osteoarthritis (“ROSA”) (CIV ID: CIV-19-01-026964)

Contura has also sponsored a multi-centre, randomised, double-blind, clinical trial to compare the effectiveness of intra-articular polyacrylamide hydrogel (Arthrosamid®) and a hyaluronic acid (HA), Synvisc-One® to safely induce symptomatic benefit in subjects with knee OA.

A total of 239 subjects were randomised into this study, which was conducted in accordance with the Declaration of Helsinki and ISO 14155 at 3 sites in Denmark.

The study was initially planned to end after 1 year but was extended to follow the participants of the Arthrosamid® arm for up to 5 years through yearly, unblinded follow-up visits.

Three-year data are available from this study. One-year data were published in 2022 and 2023 [9-11].

5.2.5.1 Clinical investigation design

Clinical investigation objectives and endpoints

The primary objective of the main study was to compare the effectiveness of one injection of 6 mL intra-articular Arthrosamid® with that of one injection of 6 mL Synvisc-One® on pain over 6 months in subjects with mild to severe knee OA.

The primary endpoint of the main study was:

- Change in the WOMAC pain subscale from baseline to week 26

The secondary objectives were to compare the effectiveness of 6 mL Arthrosamid® with that of 6 mL Synvisc-One® on knee OA-related subject reported symptoms, function during activities of daily living, and quality of life (QoL) over 1, 3, 6 and 12 months.

The secondary endpoints of the main study were change from baseline to weeks 4, 12, 26, and 52 in:

- The WOMAC pain, stiffness and physical function related subscales
- PGA
- OMERACT-OARSI response
- EuroQoL-5D-5L

The objectives of the extension study are to evaluate the long-term effectiveness of 6 mL Arthrosamid® on knee OA and to evaluate time to start of other treatments for knee osteoarthritis.

The endpoints of the extension study are:

- Change from baseline in WOMAC pain, stiffness, and physical function
- Change from baseline in subject's global assessment of impact of osteoarthritis (PGA)
- Time from treatment until target knee replacement
- Time from treatment to start of any other treatment for knee osteoarthritis

The safety objectives (for both main and extension study) are to evaluate safety of Arthrosamid® in subjects with knee OA.

Safety endpoints (for both main and extension study) are:

- Incidence of adverse events (AEs)
- Incidence of adverse device effects (ADEs)

Study analysis set definitions

The following analysis sets were used in the analysis of the data:

ITT analysis set: The ITT analysis set contains all randomised subjects irrespective of whether the subject actually received study intervention or the subject's compliance with the study protocol. The subjects were included as randomised.

Full analysis set (FAS): The FAS contains all randomised subjects who received study treatment (i.e. Arthrosamid® or Synvisc-One®) and who have WOMAC pain subscale assessed at baseline and at least one post-baseline assessment. The subjects were included as randomised.

Per protocol (PP) analysis set: The PP analysis set was defined as all subjects in the FAS, meeting all inclusion criteria, and who do not have any protocol deviations of clinical or statistical significance.

Safety analysis set: The safety analysis set is defined as subjects who have received a study treatment. The subjects will be included as actually treated.

Study treatment

All subjects were given a local anaesthetic prior to the injection, and all received prophylactic treatment with azithromycin and moxifloxacin.

Subject population

Inclusion criteria:

A patient was eligible for inclusion in the trial if all of the following criteria applied:

1. Signed and dated informed consent form.
2. Stated willingness to comply with all study procedures and availability for the duration of the study.
3. Male or female, aged ≥ 18 years.
4. Clinical diagnosis of knee OA according to American College of Rheumatology criteria.
5. Definite radiographic OA in the most symptomatic knee (target knee) at mild to severe stage (Kellgren-Lawrence 2-4).
6. If subject is using analgesics for knee OA the dose should have been stable for the past four weeks.
7. NRS (11 points (0-10) pain intensity numerical rating scale) ≥ 4 in target knee, during the past week when walking.
8. Body Mass Index (BMI) between 20-35 kg/m².
9. For females of reproductive potential: use of adequate contraception must be used throughout the trial.

Exclusion criteria:

A patient was not eligible for inclusion in the trial if any of the following criteria applied:

1. Female participants who are pregnant, lactating or planning pregnancy during the course of the clinical investigation.
2. Previous intra-articular injection of PAAG in the target knee.

3. Previous intra-articular injection with hyaluronic acid or derivatives in target knee in the previous 6 months.
4. Significant valgus/varus deformity of the knee, ligamentous laxity or meniscal instability.
5. Diseases in target knee other than OA.
6. Intra-articular injection of any substance other than hyaluronic acid (e.g. corticosteroids) in the target knee within the last 3 months.
7. Acute serious infection of any region that required hospital care or intravenous antibiotic treatment within the last 30 days, or oral antibiotic treatment within the last 14 days.
8. Skin disease or infections in the area of the injection site.
9. Infected or severely inflamed knees.
10. History of sepsis in any joint or any clinical concern for an infectious process in the target knee.
11. History of surgery in the target knee within the past 6 months.
12. Symptomatic OA of the hips, spine or ankle, that interferes with the evaluation of the target knee.
13. Planned surgery on any lower extremity.
14. Clinically significant venous or lymphatic stasis present in the legs.
15. Suffering from any unstable or severe cardio-vascular disease.
16. Any other contraindication to intra-articular injection.
17. Any foreign material in the target joint.
18. Any significant medical condition (e.g. significant psychiatric or neurological disorders or active alcohol/drug abuse) that is unstable/poorly controlled or other factors that may interfere with study participation.
19. Treatment with systemic steroids.
20. Significant change in physiotherapy in lower extremities related to OA within the previous month.
21. Fibromyalgia.
22. Inflammatory or other disease/condition which may affect the knee joint (e.g. rheumatoid arthritis, psoriatic arthritis and ankylosing spondylitis among others).
23. Haemophilia.
24. Any other condition that in the opinion of the investigator puts a potential participant at risk or otherwise precludes participation in the investigation.
25. Known allergic reactions to components of Synvisc-One (avian protein).
26. Known allergic reactions to antibiotics (azithromycin and moxifloxacin) or local anaesthesia.
27. Participation in any experimental device study within 6 months prior to enrolment, or participation in an experimental drug study within 1 month prior to enrolment.

5.2.5.2 Results

Disposition of subjects

258 subjects were screened, and 239 subjects were randomised into the study (120 to Synvisc-one and 119 to Arthrosamid®). 215 subjects completed the main study. Two subjects randomised to receive Synvisc-one were mistakenly treated with Arthrosamid®. Therefore, the safety analysis set consists of 121 subjects in the Arthrosamid® arm and 118 subjects in the Synvisc-one arm.

91 subjects entered the extension phase, from which 75 completed the 3-year visit.

Effectiveness endpoints

Primary effectiveness endpoint – main study

Descriptive results:

The scores for the transformed WOMAC pain subscale at baseline, weeks 4, 12, 26 and 52 are summarised for the FAS analysis set together with the respective changes from baseline in [Table 14](#). The mean transformed WOMAC pain score was lowest for both arms at 12 weeks. At 26 and 52 weeks the mean transformed WOMAC pain score appeared to increase slightly in the Synvisc-one arm compared to a more stable response in the Arthrosamid® arm.

The primary endpoint analysis was repeated for the ITT and PP analyses sets. Similar results were seen for these sets.

Table 14. Summary of transformed WOMAC pain subscale (0-100) – Full analysis set.

	Arthrosamid®		Synvisc-One	
	Actual values	Change from baseline	Actual values	Change from baseline
Baseline				
N	119		120	
Mean (SD)	45.1 (13.4)		46.5 (13.3)	
Median	45.0		50.0	
Min - Max	10 - 75		15 - 75	
Week 4				
N	119	119	120	120
Mean (SD)	32.6 (18.0)	-12.5 (18.8)	30.2 (15.9)	-16.3 (16.7)
Median	30.0	-10.0	30.0	-15.0
Min - Max	0 - 100	-65 - 75	0 - 65	-55 - 15
Week 12				
N	115	115	119	119
Mean (SD)	25.9 (16.7)	-19.4 (19.2)	28.6 (15.7)	-18.1 (17.2)
Median	25.0	-20.0	30.0	-20.0
Min - Max	0 - 65	-65 - 25	0 - 60	-65 - 20
Week 26				
N	114	114	117	117
Mean (SD)	26.8 (17.9)	-18.4 (19.0)	31.0 (18.4)	-15.4 (18.7)
Median	25.0	-15.0	30.0	-15.0
Min - Max	0 - 75	-65 - 20	0 - 75	-60 - 30
Week 52				
N	107	107	109	109
Mean (SD)	26.3 (19.5)	-18.2 (21.4)	31.9 (18.1)	-14.9 (17.2)
Median	25.0	-15.0	35.0	-10.0
Min - Max	0 - 75	-65 - 20	0 - 65	-65 - 20

N: Number of subjects, SD: Standard deviation.

Table 15. Analyses of change from baseline to week 26 in transformed WOMAC pain subscale (0-100) – by analysis set.

	N	LSMean (95% CI)	Treatment difference (95% CI)	Non-inferiority* (yes/no)	p-value
Full analysis set					
Synvisc-One	117	-14.83 (-17.96; -11.69)	3.59 (-0.87; 8.05)	Yes	NA
Arthroamid®	114	-18.41 (-21.58; -15.25)			
ITT analysis set					
Synvisc-One	120	-14.77 (-17.88; -11.65)	3.63 (-0.79; 8.06)	Yes	NA
Arthroamid®	119	-18.40 (-21.53; -15.27)			
PP analysis set					
Synvisc-One	99	-16.42 (-19.72; -13.12)	2.23 (-2.45; 6.91)	Yes	NA
Arthroamid®	98	-18.65 (-21.96; -15.34)			

CI: confidence interval; N: Number of subjects, LSMean: Least squares mean.

*If the lower bound of the 95% CI is > -9 the objective of non-inferiority has been met. If the lower bound of the 95% CI is > 0 superiority is declared, and the p-value for the test of superiority is presented.

Primary analysis results:

The primary analysis of the change from baseline in the transformed WOMAC pain subscale (0-100) from baseline to week 26 was performed using a mixed model for repeated measures (MMRM) including fixed, categorical effects of treatment, week, treatment-by-week interaction and site, as well as the baseline value and baseline-by-week interaction as covariates.

The primary analysis (FAS) showed a treatment difference at 26 weeks of 3.59 (95% CI: -0.87; 8.05) in favour of the Arthroamid® arm over the Synvisc-one arm (Table 15). The non-inferiority objective was met, but superiority was not reached.

For the ITT and PP analysis sets, the treatment differences were also in favour of the Arthroamid® arm over the Synvisc-one arm (Table 15) and the non-inferiority objective was met.

Secondary effectiveness endpoints – main study

WOMAC pain subscale

The analysis of the change from baseline in the transformed WOMAC pain subscale (0-100) from baseline to weeks 4, 12 and 52 was performed using the same MMRM as for the primary endpoint. The analysis showed no significant difference in the change from baseline transformed WOMAC pain subscale scores between the Arthroamid® and Synvisc-one arms (Table 16), although the difference at week 52 approached significance (p=0.0572).

Table 16. Analysis of change from baseline in transformed WOMAC pain subscale (0-100) –ITT analysis set.

	N	LSMean (95% CI)	Treatment difference (95% CI)	p-value
Week 4				
Synvisc-One	120	-15.90 (-18.80; -13.00)		
Arthrosamid®	119	-12.93 (-15.85; -10.02)	-2.97 (-7.08; 1.15)	0.1567
Week 12				
Synvisc-One	119	-17.37 (-20.20; -14.54)		
Arthrosamid®	115	-19.51 (-22.37; -16.65)	2.14 (-1.89; 6.16)	0.2965
Week 52				
Synvisc-One	109	-13.32 (-16.66; -9.98)		
Arthrosamid®	107	-17.93 (-21.29; -14.56)	4.61 (-0.14; 9.36)	0.0572

CI: confidence interval; N: Number of subjects, LSMean: Least squares mean.

WOMAC stiffness subscale

The analysis of the change from baseline in the transformed WOMAC stiffness subscale (0-100) was performed using the same MMRM as for the primary endpoint. The analysis showed a statistically significant treatment difference at week 4 in favour of Synvisc-one (p = 0.0315) while treatment differences in favour of Arthrosamid® at weeks 12, 26 and 52 were not significant ([Table 17](#)).

Table 17. Analysis of change from baseline in transformed WOMAC stiffness subscale (0-100) - ITT analysis set.

	N	LSMean (95% CI)	Treatment difference (95% CI)	p-value
Week 4				
Synvisc-One	120	-18.52 (-22.51; -14.53)		
Arthrosamid®	119	-12.31 (-16.32; -8.31)	-6.21 (-11.86; -0.55)	0.0315
Week 12				
Synvisc-One	119	-17.76 (-21.50; -14.03)		
Arthrosamid®	115	-20.14 (-23.92; -16.36)	2.38 (-2.94; 7.69)	0.3792
Week 26				
Synvisc-One	117	-13.96 (-18.13; -9.80)		
Arthrosamid®	114	-18.59 (-22.79; -14.38)	4.62 (-1.30; 10.54)	0.1254
Week 52				
Synvisc-One	109	-12.86 (-17.17; -8.56)		
Arthrosamid®	107	-17.87 (-22.21; -13.53)	5.01 (-1.11; 11.12)	0.1080

CI: confidence interval; N: Number of subjects, LSMean: Least squares mean.

WOMAC physical function subscale

The analysis of the change from baseline in the transformed WOMAC physical function subscale (0-100) was performed using the same MMRM as for the primary endpoint. The analysis showed a treatment difference in favour of Synvisc-one at week 4 and treatment differences in favour of Arthrosamid® at weeks 12, 26 and 52. None of the differences were statistically significant ([Table 18](#)).

Table 18. Analysis of change from baseline in transformed WOMAC physical function subscale (0-100) - ITT analysis set.

	N	LSMean (95% CI)	Treatment difference (95% CI)	p-value
Week 4				
Synvisc-One	120	-17.22 (-20.29; -14.15)		
Arthroamid®	119	-13.46 (-16.54; -10.37)	-3.77 (-8.12; 0.59)	0.0896
Week 12				
Synvisc-One	119	-18.12 (-20.91; -15.32)		
Arthroamid®	115	-21.36 (-24.19; -18.54)	3.25 (-0.73; 7.22)	0.1091
Week 26				
Synvisc-One	117	-18.13 (-21.30; -14.96)		
Arthroamid®	114	-18.90 (-22.10; -15.71)	0.77 (-3.73; 5.27)	0.7357
Week 52				
Synvisc-One	109	-15.21 (-18.60; -11.82)		
Arthroamid®	107	-17.74 (-21.16; -14.33)	2.54 (-2.28; 7.35)	0.3006

CI: confidence interval; N: Number of subjects, LSMean: Least squares mean.

Patient global assessment

The patient global assessment was to the question “How much does the knee osteoarthritis (target knee) as a whole affect your life at present?”. The response was indicated on a 100 mm visual analogue scale between “Not at all” (0) and “The worst possible” (100).

The analysis of the change from baseline in the patient global assessment was performed using the same MMRM as for the primary endpoint. The analysis showed a treatment difference in favour of Synvisc-one at week 4 and treatment differences in favour of Arthroamid® at weeks 12, 26 and 52. None of the differences were statistically significant (Table 19).

Table 19. Analysis of change from baseline in patient global assessment - ITT analysis set.

	N	LSMean (95% CI)	Treatment difference (95% CI)	p-value
Week 4				
Synvisc-One	119	-22.24 (-26.29; -18.20)		
PAAG-OA	118	-16.98 (-21.04; -12.92)	-5.27 (-10.99; 0.46)	0.0713
Week 12				
Synvisc-One	119	-18.42 (-22.44; -14.39)		
PAAG-OA	115	-22.50 (-26.57; -18.42)	4.08 (-1.65; 9.81)	0.1617
Week 26				
Synvisc-One	117	-14.45 (-18.91; -9.98)		
PAAG-OA	114	-16.79 (-21.30; -12.27)	2.34 (-4.01; 8.69)	0.4687
Week 52				
Synvisc-One	109	-13.51 (-18.12; -8.90)		
PAAG-OA	107	-17.53 (-22.18; -12.88)	4.02 (-2.53; 10.57)	0.2275

CI: confidence interval; N: Number of subjects, LSMean: Least squares mean.

OMERACT-OARSI response

The analysis of the proportion of positive responders to the OMERACT-OARSI criteria was performed using a repeated measures logistic model treatment. The analysis showed a treatment ratio (odds ratio) in favour of Synvisc-one at week 4 that was not statistically significant (Table 20). At weeks 12, 26 and 52, the treatment ratio was close to 1, indicating no difference between the 2 groups.

Table 20. Analysis of OMERACT-OARSI criteria - ITT analysis set.

	E/n	(%)	Synvisc-One / PAAG-OA Odds ratio (95% CI)	p-value
Week 4				
Synvisc-One	74 / 120	(61.7)		
PAAG-OA	64 / 119	(53.8)	1.38 (0.82; 2.31)	0.2200
Week 12				
Synvisc-One	77 / 119	(64.7)		
PAAG-OA	79 / 115	(68.7)	0.87 (0.50; 1.49)	0.5994
Week 26				
Synvisc-One	76 / 117	(65.0)		
PAAG-OA	76 / 114	(66.7)	0.95 (0.55; 1.62)	0.8466
Week 52				
Synvisc-One	62 / 109	(56.9)		
PAAG-OA	60 / 107	(56.1)	1.01 (0.59; 1.73)	0.9612

CI: confidence interval; E: Number of OMERACT-OARSI responders; n: number with non-missing values.

EuroQoL-5D-5L

The answers from each of the 5 dimensions from EQ-5D-5L were combined into one measurement termed the EQ-5D-5L. The analysis of the change from baseline in the EQ-5D-5L health state showed a statistically significant treatment difference in favour of Synvisc-one at week 4 while the results at weeks 12, 26 and 52 were similar between the Arthroamid® arm and the Synvisc-one arm (Table 21).

Table 21. Analysis of change from baseline in in EQ-5D-5L health state - ITT analysis set.

	N	LSMean (95% CI)	Treatment difference (95% CI)	p-value
Week 4				
Synvisc-One	120	0.08 (0.06; 0.10)		
PAAG-OA	119	0.04 (0.02; 0.07)	0.03 (0.01; 0.06)	0.0215
Week 12				
Synvisc-One	119	0.07 (0.05; 0.09)		
PAAG-OA	115	0.08 (0.07; 0.10)	-0.01 (-0.04; 0.02)	0.3930
Week 26				
Synvisc-One	117	0.06 (0.04; 0.08)		
PAAG-OA	114	0.06 (0.04; 0.08)	0.00 (-0.03; 0.04)	0.7976
Week 52				
Synvisc-One	109	0.05 (0.03; 0.08)		
PAAG-OA	107	0.06 (0.04; 0.09)	-0.01 (-0.05; 0.02)	0.5377

CI: confidence interval; N: Number of subjects, LSMean: Least squares mean.

Effectiveness endpoints – extension study

For the extension study effectiveness endpoints, the MMRM model was adjusted to examine the change from baseline for the Arthrosamid® arm. In addition, each of the WOMAC endpoints and PGA at the year 2 and 3 visit were analysed using the BOCF approach for subjects lacking a year 2 or 3 measurement. This analysis was performed using an ANCOVA model with treatment and site as fixed effects and baseline value as covariates.

WOMAC pain subscale

The analysis of the change from baseline in the transformed WOMAC pain subscale (0-100) showed a significant reduction ($p < 0.0001$) of 13.1 points from baseline to year 3 for the full ITT analysis set and 15.3 points for the year-2 extension subset (Table 22) indicating reduced pain in the treated knee. The most conservative approach, BOCF, also resulted in a significant reduction ($p < 0.0001$) of approximately 10.8 points from baseline to year 3 (Table 22).

Table 22. Analyses of change from baseline to 3 years in transformed (0-100) WOMAC pain subscale.

	Number of participants	LSMean (95% CI)	p-value
Year 2			
Planned analysis (MMRM)	87	-16.8 (-20.4; -13.2)	<.0001
Extension participants (MMRM)	87	-19.1 (-22.7; -15.5)	<.0001
Year 3			
Planned analysis (MMRM)	75	-13.1 (-17.9; -8.4)	<0.0001
Extension participants (MMRM)	75	-15.3 (-19.9; -10.6)	<0.0001
Baseline carried forward (ANCOVA)	119	-10.8 (-14.3; -7.4)	<0.0001

CI: Confidence interval; N: Number of subjects, LSMean: Least squares mean; WOMAC: Western Ontario and McMaster Universities Osteoarthritis Index.

WOMAC stiffness subscale

The analysis of the change from baseline in the transformed WOMAC stiffness subscale (0-100) showed a significant reduction ($p < 0.0001$) of approximately 14.4 points from baseline to year 3 for the full ITT analysis set and 17.4 points for the year-2 extension subset (Table 23), indicating reduced stiffness in the treated knee. The most conservative approach, BOCF, also resulted in a significant reduction ($p < 0.0001$) of approximately 12.4 points from baseline to year 3 (Table 23).

Table 23. Analyses of change from baseline to 3 years in transformed (0-100) WOMAC stiffness subscale.

	Number of participants	LSMean (95% CI)	p-value
Year 2			
Planned analysis (MMRM)	87	-16.1 (-20.8; -11.3)	<.0001
Extension participants (MMRM)	87	-18.6 (-23.5; -13.8)	<.0001
Year 3			
Planned analysis (MMRM)	75	-14.4 (-19.7; -9.1)	<0.0001
Extension participants (MMRM)	75	-17.4 (-22.7; -12.1)	<0.0001
Baseline carried forward (ANCOVA)	119	-12.4 (-16.5; -8.3)	<0.0001

CI: Confidence interval; N: Number of subjects, LSMean: Least squares mean; WOMAC: Western Ontario and McMaster Universities Osteoarthritis Index.

WOMAC physical function subscale

The analysis of the change from baseline in the transformed WOMAC physical function subscale (0-100) showed a significant reduction ($p < 0.0001$) of approximately 11.8 points from baseline to year 3 for the full ITT analysis set (Table 24) and 14.3 points for the year-2 extension subset (Table 24), indicating an improved physical function in the treated knee. The most conservative approach, BOCF, also resulted in a significant reduction ($p < 0.0001$) of approximately 10.4 points from baseline to year 3 (Table 24).

Table 24. Analyses of change from baseline to 3 years in transformed (0-100) WOMAC physical function subscale.

	Number of participants	LSMean (95% CI)	p-value
Year 2			
Planned analysis (MMRM)	87	-16.1 (-19.9; -12.4)	<.0001
Extension participants (MMRM)	87	-18.5 (-22.2; -14.8)	<.0001
Year 3			
Planned analysis (MMRM)	75	-11.8 (-16.5; -7.2)	<0.0001
Extension participants (MMRM)	75	-14.3 (-18.8; -9.8)	<0.0001
Baseline carried forward (ANCOVA)	119	-10.4 (-13.8; -7.1)	<0.0001

CI: Confidence interval; N: Number of subjects, LSMean: Least squares mean; WOMAC: Western Ontario and McMaster Universities Osteoarthritis Index.

Patient global assessment

The analysis of the change from baseline in the patient global assessment showed a significant reduction ($p < 0.0001$) of 18.3 points from baseline to year 3 for the full ITT analysis set and 20.0 points for the year-2 extension subset (Table 25). The most conservative approach, BOCF, also resulted in a significant reduction ($p < 0.0001$) of 14.0 points from baseline to year 3 (Table 25).

Table 25. Analysis of change from baseline to 3 years in patient global assessment.

	Number of participants	LSMean (95% CI)	p-value
Year 2			
Planned analysis (MMRM)	87	-20.3 (-25.0; -15.5)	<.0001
Extension participants (MMRM)	87	-22.4 (-27.3; -17.6)	<.0001
Year 3			
Planned analysis (MMRM)	75	-18.3 (-24.2; -12.4)	<0.0001
Extension participants (MMRM)	75	-20.0 (-25.8; -14.2)	<0.0001
Baseline carried forward (ANCOVA)	119	-14.0 (-18.4; -9.6)	<0.0001

CI: confidence interval; N: Number of subjects, LSmean: Least squares mean.

New treatment for target knee

Up to the 3-year visit, 32 subjects received a new treatment ([Table 26](#)). The median time to new treatment was 375.5 days.

Of the 32 subjects, 13 had a knee replacement with a median duration to knee replacement of 728 days ([Table 26](#)).

Table 26. Summary of new treatment (incl. knee replacement) for target knee OA - up to and including 3-year extension–safety analysis set.

Arthrosamid®	
Time to new treatment (incl. replacement) for knee OA (days)	
N	32
Mean (SD)	519.1 (401.0)
Median	375.5
Min - Max	1 - 1063
Time to target knee replacement (days)	
N	13
Mean (SD)	708.9 (312.1)
Median	728.0
Min - Max	212 - 1063

N: Number of subjects, SD: Standard deviation, OA: Osteoarthritis.

Safety endpoints

Adverse events

In the main study period, 55 subjects (45.5%) in the Arthrosamid® arm reported a total of 99 AEs compared to 78 AEs reported by 49 subjects (41.5%) in the Synvisc-one arm. Ten of the events were reported as SAEs, one death was reported, and one AE led to withdrawal from the study. The most frequently reported AEs were arthralgia (Arthrosamid®: 39 events in 33 subjects; Synvisc-one: 24

events in 23 subjects) and joint swelling (Arthrosamid®: 16 events in 15 subjects; Synvisc-one: 7 events in 6 subjects). The majority of the AEs were mild in intensity, although 5 events (arthritis infective, endocarditis, erysipelas, pneumonia, pulmonary embolism) in the Arthrosamid® arm and 2 events (cardiac arrest, adenocarcinoma) in the Synvisc-one arm were reported as severe. None of the SAEs were assessed as related to the study device.

Between the 1-year visit and the 3-year visit, 50 AEs were reported by 36 subjects in the Arthrosamid® arm. The most frequently reported AEs were arthralgia and joint swelling. The majority of the AEs were mild in intensity. Six of the AEs were reported as SAEs in 5 subjects. None of the 6 events (back pain, corona virus infection, death, 2 events of myocardial infarction, worsening of hip OA) were assessed as related to Arthrosamid®, and all events, except for the death, recovered without sequelae.

Adverse device effects

In the main study period, 35 subjects (28.9%) in the Arthrosamid® arm reported a total of 41 AEs compared to 13 AEs reported by 9 subjects (7.6%) in the Synvisc-one arm. The only ADEs reported in more than 1 subject were arthralgia (Arthrosamid®: 21 events in 19 subjects; Synvisc-one: 8 events in 7 subjects) and joint swelling (Arthrosamid®: 13 events in 13 subjects; Synvisc-one: 4 events in 3 subjects). The majority of the ADEs were mild in intensity, and none were reported as severe. No unanticipated serious ADEs were reported up to 12 months after treatment.

There were no ADEs reported for the extension study.

5.2.5.3 Conclusion

At 26 weeks after treatment, the effectiveness of Arthrosamid® was non-inferior to hyaluronic acid, as measured by the WOMAC pain subscale.

At 52 weeks after treatment, the effectiveness of Arthrosamid® was numerically superior to hyaluronic acid but not statistically significantly different.

The extension study provides valid scientific evidence that Arthrosamid® remains safe and appears effective for its intended use 3 years after injection.

5.3 Other clinical data

Searches for literature containing clinical data with direct relevance to the safety and performance of Arthrosamid® have been conducted covering the period until and including 2023.DEC.31. A summary of the identified publications is provided below.

Author, publication date	Title and type of study	Objective	Patients and treatments	Follow-up time	Effectiveness and safety conclusions
Christensen, L.H. et al., 2016 [12]	Histological appearance of the synovial membrane after treatment of knee osteoarthritis with polyacrylamide gel Injections: A case report Case report	Histological investigation of the synovial membrane after polyacrylamide hydrogel injection	1 knee OA patient was treated with Contura polyacrylamide hydrogel	9 months	Performance: Histological evidence of gel integration and augmentation of the thickness of the synovial membrane 9 months after injection. Areas of tissue network were seen within the integrated gel layer, either as a band immediately beneath the surface cells or as fine tissue strands inside the gel. Safety: No histological signs of necrosis, calcifications or fibrinous exudate.
Maulana, R. et al., 2022 [13]	Reduction in patellofemoral bone marrow lesions following single Arthrosamid intra-articular injection of polyacrylamide hydrogel (ipaag) in the treatment of advanced osteoarthritis Case report	To present the first case of a reduction in patellofemoral bone lesions on imaging following a single injection of Arthrosamid® in a patient with knee OA	1 knee OA patient was treated with Arthrosamid®	4 months	Performance: An 81-year-old woman with OA of the patellofemoral joint was injected with a single dose of 6 mL Arthrosamid®. After 4 months, the patient reported marked improvements in terms of greater mobility and reduced pain, and MRI demonstrated a reduction in the size of the bone marrow lesions seen at MRI prior to treatment. Safety: The injection of Arthrosamid® in this case was uneventful.

5.4 Summary of the clinical performance and safety

Arthrosamid® is an injectable, non-absorbable, non-degradable and non-migratory hydrogel intended to be used for symptomatic treatment of adult patients with knee osteoarthritis. Arthrosamid® is expected to provide relief of symptoms for 5 years, but will remain in the knee for the entire lifetime of the patient as an integral part of the synovial tissue of the inner joint capsule of the knee. Long-term clinical data on the stability of Contura's polyacrylamide hydrogel is available from a 10-year follow-up

study including 104 HIV-patients injected with an average of 6 mL of Aquamid®² for treatment of facial lipoatrophy. At follow-up (10 years), no patients presented with migration of the hydrogel, and the majority of patients were “highly satisfied” (74.8%) or “satisfied” (23.4%) with the cosmetic result [14]. Stable appearance of hydrogel was also observed in an 8-year follow-up study of 25 women who had been treated with Bulkamid® Hydrogel³ for stress urinary incontinence [15].

The clinical benefits of Arthrosamid® in patients with knee osteoarthritis has been documented through clinical investigations and market experience. Collectively, three clinical investigations (described in sections 5.2.1, 5.2.4 and 5.2.5) and one case report [13] have found short-term (from 4 weeks) and long-term (up to 4 years) improvements in WOMAC pain, stiffness and physical function sub-scores after injection with Arthrosamid®. Improvements in quality of life, as measured by PGA and EQ-5D-5L, were also found (sections 5.2.4 and 5.2.5), and integration of Arthrosamid® into the synovial membrane has been established through histopathological studies (section 5.2.3 and [12]). In addition, experience from the market indicates that Arthrosamid® performs according to its intended use.

Four clinical investigations (described in sections 5.2.2, 5.2.3, 5.2.4 and 5.2.5) and two case reports found in the literature [12,13] contribute with data supporting the safety of Arthrosamid®. Most of the AEs observed in the clinical investigations were short-lived and mild in intensity. The most frequently reported AEs were arthralgia and joint swelling. Experience from the market indicates that the subject device is safe when used as intended.

The harms associated with the residual risks of Arthrosamid® are transient symptoms not requiring professional medical or surgical intervention which patients may experience during the first few weeks after the procedure (recovery period). The harms associated with the residual risks are:

- Swelling/edema.
- Worsening of arthralgia.

The risk analysis did not identify any unacceptable risks.

Overall, the risks associated with the device are considered acceptable when weighed against the benefits to the patients.

5.5 Post-market clinical follow-up

A post-market clinical follow-up (PMCF) study with the primary objective of evaluating the long-term safety of Arthrosamid® (“LUNA”) was initiated in 2021. The study is a long-term, observational, all-comers open-label, multi-national clinical investigation to evaluate long-term safety and effectiveness. Participants were treated with a single intra-articular injection of 6 mL Arthrosamid® and will be followed for 5 years.

² Aquamid® has the same chemical composition as Arthrosamid®, but slightly different visco-elastic properties.

³ Bulkamid® Hydrogel has the same chemical composition and visco-elastic properties as Arthrosamid®.

The endpoints of the study are:

- Primary endpoint: Number of adverse events.
- Secondary endpoints:
 - Change from baseline to months 6, 12, 24, 36, 48 and 60 in:
 - The WOMAC pain, stiffness, and physical function related subscales.
 - Patient Global Assessment of impact of OA.

Recruitment was finalised in 2023. 203 adult patients were treated across 9 rheumatology and orthopaedic sites: 2 sites in Denmark, 2 in France, 2 in Germany, 2 in United Kingdom and 1 in Italy. No results are available.

Other PMCF activities include systematic screening and review of scientific literature and other sources of clinical data. Vigilance data is continuously monitored, and the periodic safety update report and clinical evaluation report are updated annually.

6. Possible therapeutic alternatives

The European Society for Clinical and Economic Aspects of Osteoporosis, Osteoarthritis and Musculoskeletal Diseases (ESCEO) recently published an updated algorithm with recommendations for the management of OA [16]. The ESCEO recommendation follows a stepwise algorithm from basic health-promoting principles such as education, weight loss if overweight, and exercise programs, to end-stage disease management in the form of surgical joint replacement. In this model, Arthrosamid® can be regarded as an alternative to intra-articular HA and corticosteroids and offered to persistently symptomatic patients in whom oral pharmacological management of knee OA has been insufficient. As surgery is an invasive procedure, this is only appropriate when all previous modalities have failed, if the patient is severely symptomatic, and there is a significant loss in quality of life [16].

6.1 Intra-articular treatments

6.1.1 Corticosteroids

Local intra-articular injection with steroids may be effective in treating symptoms of inflammation, but the effect wears off quickly (2-6 weeks) [17], though there are descriptions of persistent effect for up to six months [18,19]. Repeat injections are thus required for stable pain relief, increasing the risk of complications. Side-effects of this kind of therapy include infections, damage of soft tissue and swelling [17,19,20]. Joint infections have been described with up to 1 infection for every 10.000 steroid injections [20].

6.1.2 Hyaluronic acid

There are numerous HA-based viscosupplemental products for human use on the market. The preparations have varying average molecular weight (0.5-6.0 million Da) and therefore show great variations in visco-elastic properties, in turn leading to varying effects in terms of pain reduction and functional improvement [21-24] as well as required dosage and frequency of injections [25].

The general effectiveness of HA in reducing symptoms of OA is still debated. There are several recent meta-analyses indicating a better and longer-lasting effect of HA than that of placebo or corticosteroids, with the effect lasting up to 26 weeks [26-33]. On the other hand, a meta-analysis from 2015 found no difference between HA and placebo if open-label studies were excluded from the analysis [34]. Furthermore, recent evidence has found comparable efficacy and safety profiles of HA and placebo at 6 months post-injection [35,36] and a higher risk of later undergoing knee arthroplasty among patients treated with HA than among patients treated with corticosteroids [37].

Intra-articular HA treatment is generally regarded as safe, with a safety pattern similar to that of saline injection [38,39], although an increased risk of both non-serious [40] and serious [36] adverse events has been demonstrated. As HA is known to have a limited residence time in the knee [41], the need for repeat injections further increases the risk of adverse reactions. Side-effects include short-term-swelling, fluid retention and pain [42-46].

6.1.3 Noltrex

Another polyacrylamide gel for intra-articular treatment of knee OA, sold under the brand name Noltrex, was CE-marked in 2012, but to our knowledge no longer holds a CE-mark.

Noltrex is based on cross-linked polyacrylamide ($4.0\% \pm 1.5\%$) with silver ions added (0.0001 0.0025%). Intra-articularly injected Noltrex is claimed to cover the entire joint surface and the synovial lining with a uniform layer of lubricating gel, creating a shock-absorbing effect which helps to preserve the integrity of cartilage/stop further destruction of the cartilage. The joint space is claimed to widen due to viscosupplementation, leading to mechanical protection of articular cartilages. The cross-linked polymer is claimed to have a long-lasting effect compared to HA-based products. However, as opposed to Arthrosamid®, Noltrex is eliminated from the body, and the injections are to be repeated every 6 months - 1.5 years.

6.2 Surgery

The surgical techniques that may be employed include total joint replacement, partial knee replacement, or osteotomy around the knee. Although TKA is highly successful and cost-effective [16,47], it often has adverse outcomes. Unicompartamental knee replacement (UKR) has a lower occurrence of complications compared with TKA and may be associated with a higher activity level, and better functional recovery [48]. UKR may lead to a higher revision/reoperation rate than TKA, but literature is not consistent on this finding [48-50].

Most patients have a good clinical outcome after knee replacement, but multiple studies have reported that 20% or more of patients are dissatisfied with their outcome, usually characterised by ongoing pain and poor function [47,51]. It has been reported that a 65-year-old patient will have a 7% lifetime risk of requiring revision surgery, but this risk increases substantially with treatment in younger age groups. In fact, recent evidence suggests that in the age group of those younger than 60 years the lifetime risk of revision is approximately 35% for men and 20% for women, with half the rate of revision occurring within the first 5 years of implantation [47]. Post-surgical pain, instability and stiffness, together with infection and implant loosening, have been stated to account for 80–90% of all revisions [52].

Serious complications to TKA are relatively frequent. A 30-day major complication rate of approximately 2% has been reported, with pulmonary embolism being the most prevalent major systemic complication (noted in 0.78% of the patients), and deep wound infection being the most prevalent major local complication (0.3% of the patients) [53]. A general incidence of prosthesis infection around 2% has been observed [54]. Other frequent adverse events include systemic complications such as urinary tract infection and deep venous thrombosis, and local complications such as superficial wound infection [53].

Surgical vs. non-surgical intervention (physiotherapy, dietary advice, use of insoles and pain medication) in patients with moderate to severe knee OA was investigated in a randomised clinical trial including 100 patients, published in 2015 [55]. Patients either underwent TKA followed by 12 weeks of non-surgical treatment or received only the 12 weeks of non-surgical treatment. Treatment with TKA followed by non-surgical treatment resulted in greater pain relief and functional improvement after 12 months than did non-surgical treatment alone. However, TKA was associated with a higher number of serious adverse events than was non-surgical treatment. The 2-year follow-up study concluded that TKA followed by supervised non-surgical treatment resulted in twice the improvement in pain and function after 2 years compared with non-surgical treatment with the option of TKA later [56].

7. Suggested profile and training for users

Arthrosamid® is intended to be administered by qualified physicians familiar with intra-articular injection procedures.

8. Harmonised standards and applied Common Specifications

A list of standards, including harmonised standards, applied to the design, manufacture and packaging of Arthrosamid® is found in the table below.

There are no applicable Common Specifications for the conformity assessment of Arthrosamid®.

Origin	No.	Year	Subject	Title
DS/EN ISO	10993-1	2020	Biological evaluation	Biological evaluation of medical devices – Part 1 Evaluation and testing
DS/EN ISO	10993-3	2014	Biological evaluation	Biological evaluation of medical devices – Part 3: Test for genotoxicity, carcinogenicity and reproductive toxicity
DS/EN ISO	10993-5	2009	Biological evaluation	Biological evaluation of medical devices – Part 5: Tests for in vitro cytotoxicity
DS/EN ISO	10993-6	2016	Biological evaluation	Biological evaluation of medical devices – Part 6: Tests for local effects after implantation
DS/EN ISO	10993-9	2021	Biological evaluation	Biological evaluation of medical devices – Part 9 Framework for identification and degradation of potential degradation products
DS/EN ISO	10993-10	2023	Biological evaluation	Biological evaluation of medical devices – Part 10: Tests for irritation and delayed-type hypersensitivity
DS/EN ISO	10993-11	2018	Biological evaluation	Biological evaluation of medical devices – Part 11: Tests for systemic toxicity
DS/EN ISO	10993-12	2021	Biological evaluation	Biological evaluation of medical devices - Part 12: Sample preparation and reference materials
DS/EN ISO	10993-13	2010	Biological evaluation	Biological evaluation of medical devices – Part 13: Identification and quantification of degradation products from polymeric medicinal devices

Origin	No.	Year	Subject	Title
DS/EN ISO	10993-18:2020 + A1:2023	2020/2023	Biological evaluation	Biological evaluation of medical devices – Part 18: Chemical characterisation of materials
DS/ISO/TS	10993-19	2020	Biological evaluation	Biological evaluation of medical devices – Part 19: Physico-chemical, morphological and topographical characterisation of materials
DS/EN/TR	10993-33	2015	Biological evaluation	Biological evaluation of medical devices - Part 33: Guidance on tests to evaluate genotoxicity - Supplement to ISO 10993-3
DS/EN/TR	21726	2019	Biological evaluation	Biological evaluation of medical devices - Application of the threshold of toxicological concern (TTC) for assessing biocompatibility of medical devices constituents
DS/EN	12469	2000	Cleanrooms	Biotechnology - Performance criteria for microbiological safety cabinets
EN ISO	14644-1	2015	Cleanrooms	Cleanrooms and associated controlled environments – Part 1: Classification of air cleanliness
EN ISO	14644-2	2015	Cleanrooms	Cleanrooms and Associated Controlled Environments - Part 2: Specification for Testing and Monitoring to Prove Continued Compliance with ISO 14644-1
EN ISO	14644-3:2019 + Cor. 2020-06	2019/2020	Cleanrooms	Cleanrooms and Associated Controlled Environments - Part 3: Test methods
EN ISO	14644-4	2022	Cleanrooms	Cleanrooms and Associated Controlled Environments - Part 4: Design, Construction and Start-up
EN ISO	14644-5	2004	Cleanrooms	Cleanrooms and Associated Controlled Environments - Part 5: Drift
EN ISO	14644-7	2004	Cleanrooms	Cleanrooms and Associated Controlled Environments - Part 7: Separative devices (clean air hood, gloveboxes, isolators and mini-environments)
EN ISO	14644-8	2022	Cleanrooms	Cleanrooms and Associated Controlled Environments - Part 8: Classification of airborne molecular contamination
DS/EN	17141	2020	Cleanrooms	Cleanrooms and associated controlled environments - Biocontamination control
DS/EN ISO	14155	2020	Clinical evaluation	Clinical investigation of medical devices for human subjects – Good clinical practice
EN ISO	14630	2013	Development	Non-active surgical implants – General requirements
EN ISO	14971:2019 + A11:2021	2019/2021	Development	Medical devices – Application of risk management to medical devices
DS/EN	62366-1:2015 +AC:2015 +AC:2016 +A1:2020	2015/2016/2020	Development	Medical Devices - Part 1: Application of usability engineering to medical devices
DS/ISO	639	2023	Labelling; Nomenclature; Symbols	Code for individual languages and language groups
DS/ISO IEC	7810	2019	Labelling; Nomenclature; Symbols	Identification cards - Physical characteristics
DS/EN ISO	15223-1	2021	Labelling; Nomenclature; Symbols	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements
DS/EN ISO	15223-1 Bil.1	2021	Labelling; Nomenclature; Symbols	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements - Fact Sheet 1: Graphical symbols in electronic format
DS/ISO	15223-2	2011	Labelling; Nomenclature; Symbols	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied - Part 2: Symbol Development, selection and validation
DS/EN ISO	20417	2021	Labelling; Nomenclature; Symbols	Information supplied by the manufacturer with medical devices
USP	61	2023	Manufacturing and validation	Microbiological Examination of Non-Sterile Products: microbial enumeration tests labs and test for specified organisms
USP	85	2023	Manufacturing and validation	Bacterial Endotoxins Test
USP	87	2023	Manufacturing and validation	Biological Reactivity Tests, In vitro

Origin	No.	Year	Subject	Title
USP	161	2023	Manufacturing and validation	Bacterial Endotoxin and Pyrogen Tests
USP	755	2023	Manufacturing and validation	Minimum Fill
ISO	8573-1	2010	Manufacturing and validation	Compressed Air
EN ISO/IEC	17025	2017	Manufacturing and validation	General requirements for the competence of testing and calibration laboratories
EN/ISO	80369-1	2018	Manufacturing and validation	Small-bore connectors for liquids and gases in healthcare applications - part 1: general requirements
DS/EN ISO	80369-7	2021	Manufacturing and validation	Small-bore connectors for liquids and gases in healthcare applications - Part 7: Connectors for intravascular or hypodermic applications
DS/EN ISO	80369-20	2015	Manufacturing and validation	Small-bore connectors for liquids and gases in healthcare applications - Part 20: Common test methods
Ph.Eur	0169	2024	Manufacturing and validation	Nitrates
Ph.Eur	2.1.6	2024	Manufacturing and validation	Gas Detector Tubes
Ph.Eur	2.2.11	2024	Manufacturing and validation	Distillation range
Ph.Eur	2.2.12	2024	Manufacturing and validation	Boiling point
Ph.Eur	2.2.24	2024	Manufacturing and validation	Absorption Spectrophotometry, Infrared
Ph.Eur	2.2.3	2024	Manufacturing and validation	Potentiometric Determination of pH
Ph.Eur	2.2.32	2024	Manufacturing and validation	Loss on Drying
Ph.Eur	2.2.44	2024	Manufacturing and validation	Total Organic Carbon (TOC)
Ph.Eur	2.2.6	2024	Manufacturing and validation	Refractive Index
Ph.Eur	2.4.8	2024	Manufacturing and validation	Heavy metals
Ph.Eur	2.6.12	2024	Manufacturing and validation	Microbiological Examination of Non-Sterile Products: microbial enumeration tests
Ph.eur	1238	2024	Manufacturing and validation	Air Medicinal
DS/ISO	2859-1:1999 + Amd.1:2011	1999/2011	Manufacturing and validation	Sampling procedures for inspection by attributes • Part 1: Sampling schemes indexed by acceptable quality limit (AQL) for lot-by-lot inspection
ISO	4180	2019	Packaging	Packaging- Complete, filled transport packages - General rules for the compilation of performance test schedules
DS/EN ISO	11607-1:2020 + A11:2022 + A1:2023	2020/2022/2023	Packaging	Packaging for terminally sterilised medical devices – Part 1: Requirements for materials, sterile barrier systems and packaging systems
EN ISO	11607-2:2020 + A11:2022 + A1:2023	2020/2022/2023	Packaging	Packaging for terminally sterilised medical devices – Part 2: Validation requirements for forming, sealing and assembly processes
ANSI/AAMI	ST67	2019	Sterilization	Sterilization of health care products - requirements and guidance for selecting a sterility assurance level (SAL) for products labeled "sterile"
ASTM	F 1929-15	2015	Packaging	Standard Test Method for Detecting Seal Leaks in Porous Medical Packaging by Dye Penetration
ASTM	F 1980-21	2021	Packaging	Standard Guide for accelerated Aging of Sterile Barrier Systems and Medical Devices
ASTM	F 2096 – 11 (2019)	2019	Packaging	Standard Test Method for Detecting Gross Leaks in Medical Packaging by Internal Pressurisation (Bubble Test)
ASTM	F 1886/F1886M-16	2016	Packaging	Standard Test Method for Determining Integrity of Seals for Flexible Packaging by Visual Inspection

Origin	No.	Year	Subject	Title
ASTM	F88/F 88M-23	2023	Packaging	Standard Test Method for Seal Strength of Flexible Barrier Materials
DS/EN ISO	13485:2016 + AC:2018 + A11:2021	2016/2018/2021	QA-general	Medical devices – Quality management systems – Requirements for regulatory purposes
ISO	28594	2017	Sampling procedures	Combined accept-zero sampling systems and process control procedures for product acceptance
ISO	2859-1:1999 + Cor.1::2001 + Amd 1:2011	1999/2001/2011	Sampling procedures	Sampling procedures for inspection by attributes - Part 1: Sampling schemes indexed by acceptable quality limit (AQL) for lot-by-lot inspection
DS/EN/ ISO	11138-1	2017	Sterilisation	Sterilisation of health care products - biological indicators - Part 1
DS/EN ISO	11138-3	2017	Sterilisation	Sterilisation of health care products - biological indicators - Part 3 Biological indicators for moist heat sterilisation indicators
DS/EN ISO	11138-7	2019	Sterilization	Sterilization of health care products - biological indicators – Part 7: Guidance for selection, use and interpretation of results
EN ISO	17665-1	2006	Sterilisation	Sterilisation of health care products – Moist heat – Part 1: Requirements for the development, validation and routine control of a sterilisation process for medical devices
EN	556-1:2001 + AC:2006	2001/2006	Sterilisation	Sterilisation of medical devices to be designated “STERILE” – Part 1: Requirements for terminally sterilised medical devices

9. Revision history

SSCP revision number	Date	Description of changes	Revision validated by notified body
018.SSCP.00023.03 ver. 001	2024-06-04	New document	☒ Yes Validation language: English ☐ No

A summary of the safety and clinical performance of the device, intended for patients, is given below.

10. Bibliography

1. Henriksen, M., Overgaard, A., Hartkopp, A., Bliddal, H., 2018. Intra-articular 2.5% polyacrylamide hydrogel for the treatment of knee osteoarthritis: an observational proof-of-concept cohort study. *Clinical and experimental rheumatology* 36, 1082-1085.
2. Overgaard, A., Bliddal, H., Henriksen, M., 2019. Safety of Intra-Articular Polyacrylamide Hydrogel for the Treatment of Knee Osteoarthritis Symptoms: A Retrospective Case Series. *Clin Ortho Adv Res J: COARJ-100001*.
3. Christensen, L., 2017. Case series: Report on a prospective histopathological study on biopsies from seven total arthroplasty knees (TKA) treated with polyacrylamide hydrogel (PAAG). Internal document.
4. Bliddal, H., Overgaard, A., Hartkopp, A., Beier, J., Conaghan, P.G., Henriksen, M., 2021. Polyacrylamide Hydrogel Injection for Knee Osteoarthritis: A 6 Months Prospective Study. *J Orthop Res Ther* 6, 1188.
5. Bliddal, H., Overgaard, A., Hartkopp, A., Beier, J., Conaghan, P.G., Henriksen, M., 2021. Polyacrylamide Hydrogel Injection for Knee Osteoarthritis: Results of a 52 Week Prospective Study. *Osteoarthritis and Cartilage* 29, S278.
6. Bliddal, H., Beier, J., Hartkopp, A., Conaghan, P.G., Henriksen, M., 2022. A Prospective Study of Polyacrylamide Hydrogel for Knee Osteoarthritis: Results From 2 Years After Treatment. *Osteoarthritis and Cartilage* 30, S371-S372.
7. Henriksen, M., Beier, J., Hartkopp, A., Conaghan, P.G., Bliddal, H., 2023. 3 year results from a prospective study of polyacrylamide hydrogel for knee osteoarthritis. *Osteoarthritis and Cartilage* 31, P682-3.
8. Bliddal, H., Beier, J., Hartkopp, A., Conaghan, P.G., Henriksen, M., 2023. AB1185 Polyacrylamide hydrogel for the treatment of knee osteoarthritis: 3 year follow up results of a prospective clinical study. *Annals of the Rheumatic Diseases* 82:1822-1823.
9. Bliddal, H., Beier, J., Hartkopp, A., Conaghan, P.G., Henriksen, M., 2022. One-Year Performance of Polyacrylamide Hydrogel vs. Hyaluronic Acid: A Randomised Controlled Study. *Osteoarthritis and Cartilage* 30, S370-S371.
10. Bliddal, H., Beier, J., Hartkopp, A., Conaghan, P.G., Henriksen, M., 2022. One-Year Performance of Polyacrylamide Hydrogel vs. Hyaluronic Acid in Age, BMI and Kellgren-Lawrence Subgroups: a Subgroup Analysis of a Randomised Study. *Osteoarthritis and Cartilage* 30, S373-S374.
11. Bliddal, H., Beier, J., Hartkopp, A., Conaghan, P.G., Henriksen, M., 2023. Polyacrylamide hydrogel versus hyaluronic acid in knee osteoarthritis: a subgroup analysis of a randomized study. *Orthopaedic Proceedings* 105-B(SUPP 13):81.

12. Christensen, L., Daugaard, S., 2016. Histological Appearance of the Synovial Membrane after Treatment of Knee Osteoarthritis with Polyacrylamide Gel Injections: A Case Report. *Journal of Arthritis* 5, 217.
13. Maulana, R., Cole, A., Lee, P.Y.F., 2022. Reduction in Patellofemoral Bone Marrow Lesions Following Single Arthrosamid Intra -Articular Injection of Polyacrylamide Hydrogel (ipaag) in the Treatment of Advanced Osteoarthritis. *J Arthritis* 11, 24-26.
14. Negredo, E., Puig, J., Ornelas, A., Echeverria, P., Bonjoch, A., Estany, C., Higuera, C., Gonzalez-Mestre, V., Clotet, B., 2015. Ten-Year Safety with Polyacrylamide Gel Used to Correct Facial Lipoatrophy in HIV-Infected Patients. *AIDS Res Hum Retroviruses* 31, 817-821.
15. Mouritsen, L., Lose, G., Moller-Bek, K., 2014. Long-term follow-up after urethral injection with polyacrylamide hydrogel for female stress incontinence. *Acta Obstet Gynecol Scand* 93, 209-212.
16. Bruyere, O., Honvo, G., Veronese, N., Arden, N.K., Branco, J., Curtis, E.M., Al-Daghri, N.M., Herrero-Beaumont, G., Martel-Pelletier, J., Pelletier, J.P., Rannou, F., Rizzoli, R., Roth, R., Uebelhart, D., Cooper, C., Reginster, J.Y., 2019. An updated algorithm recommendation for the management of knee osteoarthritis from the European Society for Clinical and Economic Aspects of Osteoporosis, Osteoarthritis and Musculoskeletal Diseases (ESCEO). *Seminars in arthritis and rheumatism* 49, 337-350.
17. Najm, A., Alunno, A., Gwinnutt, J.M., Weill, C., Berenbaum, F., 2021. Efficacy of intra-articular corticosteroid injections in knee osteoarthritis: A systematic review and meta-analysis of randomized controlled trials. *Joint Bone Spine* 88, 105198.
18. Tammachote, N., Kanitnate, S., Yakumpor, T., Panichkul, P., 2016. Intra-Articular, Single-Shot Hylan G-F 20 Hyaluronic Acid Injection Compared with Corticosteroid in Knee Osteoarthritis: A Double-Blind, Randomized Controlled Trial. *The Journal of bone and joint surgery. American* volume 98, 885-892.
19. Wang, C.P., Lee, W.C., Hsieh, R.L., 2022. Effects of Repeated Co-Injections of Corticosteroids and Hyaluronic Acid on Knee Osteoarthritis: A Prospective, Double-Blind Randomized Controlled Trial. *Am J Med.* 135, 641-649.
20. Marsland, D., Mumith, A., Barlow, I.W., 2014. Systematic review: the safety of intra-articular corticosteroid injection prior to total knee arthroplasty. *The Knee* 21, 6-11.
21. Gomis, A., Pawlak, M., Balazs, E.A., Schmidt, R.F., Belmonte, C., 2004. Effects of different molecular weight elastoviscous hyaluronan solutions on articular nociceptive afferents. *Arthritis and rheumatism* 50, 314-326.
22. Altman, R.D., Bedi, A., Karlsson, J., Sancheti, P., Schemitsch, E., 2016. Product Differences in Intra-articular Hyaluronic Acids for Osteoarthritis of the Knee. *The American journal of sports medicine* 44, 2158-2165.

23. Sun, S.F., Hsu, C.W., Lin, H.S., Liou, I.H., Chen, Y.H., Hung, C.L., 2017. Comparison of Single Intra-Articular Injection of Novel Hyaluronan (HYA-JOINT Plus) with Synvisc-One for Knee Osteoarthritis: A Randomized, Controlled, Double-Blind Trial of Efficacy and Safety. *The Journal of bone and joint surgery. American volume* 99, 462-471.
24. Zhao, H., Liu, H., Liang, X., Li, Y., Wang, J., Liu, C., 2016. Hylan G-F 20 Versus Low Molecular Weight Hyaluronic Acids for Knee Osteoarthritis: A Meta-Analysis. *BioDrugs : clinical immunotherapeutics, biopharmaceuticals and gene therapy* 30, 387-396.
25. Balazs, E.A., 2004. Viscosupplementation for treatment of osteoarthritis: from initial discovery to current status and results. *Surgical technology international* 12, 278-289.
26. He, W.W., Kuang, M.J., Zhao, J., Sun, L., Lu, B., Wang, Y., Ma, J.X., Ma, X.L., 2017. Efficacy and safety of intraarticular hyaluronic acid and corticosteroid for knee osteoarthritis: A meta-analysis. *International journal of surgery (London, England)* 39, 95-103.
27. Leighton, R., Akermark, C., Therrien, R., Richardson, J.B., Andersson, M., Todman, M.G., Arden, N.K., 2014. NASHA hyaluronic acid vs. methylprednisolone for knee osteoarthritis: a prospective, multi-centre, randomized, non-inferiority trial. *Osteoarthritis and cartilage / OARS, Osteoarthritis Research Society* 22, 17-25.
28. Phillips, M., Vannabouathong, C., Devji, T., Patel, R., Gomes, Z., Patel, A., Dixon, M., Bhandari, M., 2020. Differentiating factors of intra-articular injectables have a meaningful impact on knee osteoarthritis outcomes: a network meta-analysis. *Knee surgery, sports traumatology, arthroscopy: official journal of the ESSKA* 28, 3031-3039.
29. Leighton, R., Fitzpatrick, J., Smith, H., Crandall, D., Flannery, C.R., Conrozier, T., 2018. Systematic clinical evidence review of NASHA (Durolane hyaluronic acid) for the treatment of knee osteoarthritis. *Open Access Rheumatol* 10, 43-54.
30. Maheu, E., Bannuru, R.R., Herrero-Beaumont, G., Allali, F., Bard, H., Migliore, A., 2019. Why we should definitely include intra-articular hyaluronic acid as a therapeutic option in the management of knee osteoarthritis: Results of an extensive critical literature review. *Seminars in arthritis and rheumatism* 48, 563-572.
31. Maheu, E., Rannou, F., Reginster, J.Y., 2016. Efficacy and safety of hyaluronic acid in the management of osteoarthritis: Evidence from real-life setting trials and surveys. *Seminars in arthritis and rheumatism* 45, S28-33.
32. Tiendrebeogo, E., Choueiri, M., Chevalier, X., Conrozier, T., Eymard, F., 2021. Does the Presence of Neuropathic Pain Influence the Response to Hyaluronic Acid in Patients with Knee Osteoarthritis? *Cartilage* 13, 1548S-1556S.
33. Singh, H., Knapik, D.M., Polce, E.M., Eikani, C.K., Bjornstad, A.H., Gursoy, S., Perry, A.K., Westrick, J.C., Yanke, A.B., Verma, N.N., Cole, B.J., Chahla, J.A., 2022. Relative Efficacy of Intra-articular Injections in the Treatment of Knee Osteoarthritis: A Systematic Review and Network Meta-analysis. *Am J Sports Med* 50, 3140-3148.

34. Jevsevar, D., Donnelly, P., Brown, G.A., Cummins, D.S., 2015. Viscosupplementation for Osteoarthritis of the Knee: A Systematic Review of the Evidence. *The Journal of bone and joint surgery. American volume* 97, 2047-2060.
35. Ke, Y., Jiang, W., Xu, Y., Chen, Y., Zhang, Q., Xue, Q., Lin, J., Ngai, W., Nian, G., Fazeli, M.S., Xie, Y., Zhu, Z., 2021. Efficacy and safety of a single intra-articular injection of 6 ml Hylan G-F 20 compared to placebo in Chinese patients with symptomatic knee osteoarthritis : C-SOUND study, a 26-week multicenter double-blind randomized placebo-controlled trial in China. *BMC Musculoskelet Disord* 22, 428.
36. Pereira, T.V., Jüni, P., Saadat, P., Xing, D., Yao, L., Bobos, P., Agarwal, A., Hincapié, C.A., da Costa, B.R., 2022. Viscosupplementation for knee osteoarthritis: systematic review and meta-analysis. *BMJ* 378, e069722.
37. Cui, B., Chen, Y., Tian, Y., Liu, H., Huang, Y., Yin, G., Xie, Q., 2022. Effects of medications on incidence and risk of knee and hip joint replacement in patients with osteoarthritis: a systematic review and meta-analysis. *Adv Rheumatol* 62, 22.
38. Miller, L.E., Block, J.E., 2013. US-Approved Intra-Articular Hyaluronic Acid Injections are Safe and Effective in Patients with Knee Osteoarthritis: Systematic Review and Meta-Analysis of Randomized, Saline-Controlled Trials. *Clinical medicine insights. Arthritis and musculoskeletal disorders* 6, 57-63.
39. Petterson, S.C., Plancher, K.D., 2019. Single intra-articular injection of lightly cross-linked hyaluronic acid reduces knee pain in symptomatic knee osteoarthritis: a multicenter, double-blind, randomized, placebo-controlled trial. *Knee surgery, sports traumatology, arthroscopy : official journal of the ESSKA* 27, 1992-2002.
40. Miller, L.E., Bhattacharyya, S., Parrish, W.R., Fredericson, M., Bisson, B., Altman, R.D., 2021. Safety of Intra-Articular Hyaluronic Acid for Knee Osteoarthritis: Systematic Review and Meta-Analysis of Randomized Trials Involving More than 8,000 Patients. *Cartilage* 13, 351S-363S.
41. Bagga, H., Burkhardt, D., Sambrook, P., March, L., 2006. Longterm effects of intraarticular hyaluronan on synovial fluid in osteoarthritis of the knee. *J Rheumatol* 33, 946-950.
42. Xiao, J., Hu, Y., Huang, L., Huang, Z.F., Jiang, W.Z., Luo, Y.Q., Jia, M.Y., Chen, D., Shi, Z.J., 2021. Injection route affects intra-articular hyaluronic acid distribution and clinical outcome in viscosupplementation treatment for knee osteoarthritis: a combined cadaver study and randomized clinical trial. *Drug Deliv Transl Res* 11, 279-291.
43. Benazzo, F., Perticarini, L., Padolino, A., Castelli, A., Gifuni, P., Lovato, M., Manzini, C., Giordan, N., 2016. A multi-centre, open label, long-term follow-up study to evaluate the benefits of a new viscoelastic hydrogel (Hymovis(R)) in the treatment of knee osteoarthritis. *European review for medical and pharmacological sciences* 20, 959-968.
44. Bisicchia, S., Bernardi, G., Tudisco, C., 2016. HYADD 4 versus methylprednisolone acetate in symptomatic knee osteoarthritis: a single-centre single blind prospective randomised

- controlled clinical study with 1-year follow-up. *Clinical and experimental rheumatology* 34, 857-863.
45. Hermans, J., Bierma-Zeinstra, S.M.A., Bos, P.K., Niesten, D.D., Verhaar, J.A.N., Reijman, M., 2019. The effectiveness of high molecular weight hyaluronic acid for knee osteoarthritis in patients in the working age: a randomised controlled trial. *BMC musculoskeletal disorders* 20, 196.
 46. Yan, C.H., Chan, W.L., Yuen, W.H., Yung, P.S., Ip, K.Y., Fan, J.C., Chiu, K.Y., 2015. Efficacy and safety of hylan G-F 20 injection in treatment of knee osteoarthritis in Chinese patients: results of a prospective, multicentre, longitudinal study. *Hong Kong medical journal = Xianggang yi xue za zhi* 21, 327-332.
 47. Price, A.J., Alvand, A., Troelsen, A., Katz, J.N., Hooper, G., Gray, A., Carr, A., Beard, D., 2018. Knee replacement. *Lancet* 392, 1672-1682.
 48. Peng G, Liu M, Guan Z, Hou Y, Liu Q, Sun X, Zhu X, Feng W, Zeng J, Zhong Z, Zeng Y. Patellofemoral arthroplasty versus total knee arthroplasty for isolated patellofemoral osteoarthritis: a systematic review and meta-analysis. *J Orthop Surg Res.* 2021 Apr 15;16(1):264.
 49. Liddle, A.D., Judge, A., Pandit, H., Murray, D.W., 2014. Adverse outcomes after total and unicompartmental knee replacement in 101,330 matched patients: a study of data from the National Joint Registry for England and Wales. *Lancet* 384, 1437-1445.
 50. Li C, Li Z, Shi L, Gao F, Sun W. The short-term effectiveness and safety of second-generation patellofemoral arthroplasty and total knee arthroplasty on isolated patellofemoral osteoarthritis: a systematic review and meta-analysis. *J Orthop Surg Res.* 2021 Jun 2;16(1):358.
 51. Olsen, U., Lindberg, M.F., Rose, C., Denison, E., Gay, C., Aamodt, A., Brox, J.I., Skare, Ø., Furnes, O., Lee, K., Lerdal, A., 2022. Factors Correlated With Physical Function 1 Year After Total Knee Arthroplasty in Patients With Knee Osteoarthritis: A Systematic Review and Meta-analysis. *JAMA Netw Open* 5, e2219636.
 52. Carr, A.J., Robertsson, O., Graves, S., Price, A.J., Arden, N.K., Judge, A., Beard, D.J., 2012. Knee replacement. *Lancet* 379, 1331-1340.
 53. Belmont, P.J., Jr., Goodman, G.P., Waterman, B.R., Bader, J.O., Schoenfeld, A.J., 2014. Thirty-day postoperative complications and mortality following total knee arthroplasty: incidence and risk factors among a national sample of 15,321 patients. *The Journal of bone and joint surgery. American volume* 96, 20-26.
 54. Wilson, M.G., Kelley, K., Thornhill, T.S., 1990. Infection as a complication of total knee-replacement arthroplasty. Risk factors and treatment in sixty-seven cases. *The Journal of bone and joint surgery. American volume* 72, 878-883.

55. Skou, S.T., Roos, E.M., Laursen, M.B., Rathleff, M.S., Arendt-Nielsen, L., Simonsen, O., Rasmussen, S., 2015. A Randomized, Controlled Trial of Total Knee Replacement. *The New England journal of medicine* 373, 1597-1606.
56. Skou, S.T., Roos, E.M., Laursen, M.B., Rathleff, M.S., Arendt-Nielsen, L., Rasmussen, S., Simonsen, O., 2018. Total knee replacement and non-surgical treatment of knee osteoarthritis: 2-year outcome from two parallel randomized controlled trials. *Osteoarthritis and cartilage / OARS, Osteoarthritis Research Society* 26, 1170-1180.

Summary of Safety and Clinical Performance - Part II – Patients

Document revision: 001
Date issued: 2024-06-04

This Summary of Safety and Clinical Performance (SSCP) is intended to provide public access to an updated summary of the main aspects of the safety and clinical performance of Arthrosamid®. The information presented below is intended for patients or lay persons. A more extensive summary is found in the first part of this document.

The SSCP is not intended to give general advice on the treatment of a medical condition. Please contact your healthcare professional in case you have questions about your medical condition or about the use of Arthrosamid® in your specific situation. This SSCP is not intended to replace an implant card or the instructions for use to provide information on the safe use of Arthrosamid®.

1. Device identification and general information

1.1 Device trade name

Arthrosamid®, item no. 50020 (single pack) and 50021 (six-pack).

1.2 Manufacturer

Contura International A/S
Sydmarken 23
DK-2860 Soeborg
Denmark

1.3 Basic Unique Device Identification – Device Identifier (UDI-DI)

5704101Arthrosamid3K

The Basic UDI-DI is the primary identifier of a device model.

1.4 Year when the device was first CE-marked

Arthrosamid® was CE-marked in 2021.

2. Intended use of the device

2.1 Intended purpose, indication and intended patient group

Arthrosamid® is intended to be used for the treatment of knee osteoarthritis in patients, above the age of 18 years, that have failed initial attempts (e.g., exercise or weight loss programs) to alleviate symptoms.

2.2 Contraindications

Arthrosamid® should not be injected:

- If an active skin disease or infection is present at or near the injection site.
- If the joint is infected or severely inflamed.
- If the patient has previously received treatment with an absorbable intra-articular injectable (such as hyaluronic acid) that is not yet expected to be absorbed, according to manufacturer's information for the specific product.
- If the patient has previously received treatment with a different non-absorbable injectable/implant in the knee.
- If the patient has received knee replacement surgery or has any foreign material in the knee.
- If the patient has undergone knee scope surgery within the last 6 months.
- In patients with uncontrolled haemophilia (bleeder's disease) or in patients with uncontrolled anti-coagulation treatment.

A different non-absorbable implant should not subsequently be injected in the Arthrosamid®-injected knee.

3. Device description

3.1 Device description and materials in contact with patient tissues

Arthrosamid® is a hydrogel for injection into the knee cavity. It consists of 97.5% sterile water and 2.5% polyacrylamide. The hydrogel is supplied in sterile pre-filled syringes. Recommended dosage is 6 mL.


Arthrosamid® is an injectable, non-absorbable, non-degradable and non-migratory hydrogel intended to be used for symptomatic treatment of adult patients with knee osteoarthritis. Arthrosamid® is expected to provide relief of symptoms for 5 years, but will remain in the knee for the entire lifetime of the patient.

3.2 Information about medicinal substance, if any

Arthrosamid® *does not* contain any medicinal substance.

3.3 Description of how the device is achieving its intended mode of action

Arthrosamid® integrates into the synovial tissue of the inner joint capsule of the knee and thickens the synovial membrane. Arthrosamid® provides durable cushioning of the inner capsular tissue, thereby reducing pain, decreasing joint stiffness, and improving function of the knee affected by osteoarthritis. Please see below figure.

3.4 Other devices or products intended to be used in combination with the device

Antibiotics are administered prior to injection with Arthrosamid®.

Injection of Arthrosamid® is intended to be performed using a sterile needle and with the help of ultrasound.

4. Risks and warnings

Contact your healthcare professional if you believe that you are experiencing side-effects related to Arthrosamid® or its use, or if you are concerned about risks. This document is not intended to replace a consultation with your healthcare professional if needed.

Before the procedure, your physician will tell you about the intended use, expected results, contraindications, precautions, warnings and potential undesirable effects.

4.1 How potential risks have been controlled or managed

As with any invasive joint procedure, there is a risk of infection when injecting Arthrosamid®, and therefore antibiotics must be administered prior to injection.

4.2 Remaining risks and undesirable effects

Mild to moderate pain and/or mild swelling may occur during the first few weeks to months following injection. Formation or swelling of Baker's cysts (accumulation of fluid in the back of the knee) may also occur.

Infections at the injection site may occur in rare cases.

4.3 Warnings or precautions

Avoid strenuous activity (e.g., tennis, jogging or long walks) during the first weeks after the injection.

If you should need surgery or major dental work after your treatment with Arthrosamid[®], antibiotic treatment to reduce risk of infection should be considered by the treating physician.

After your treatment, inform health care professionals of your Arthrosamid[®] implant for accurate future medical assessment.

Arthrosamid[®] should be used with caution in patients:

- with acute or chronic infections in other parts of the body.
- with autoimmune disorders.
- in immunosuppressive therapy.
- with a compromised immune system.
- suffering from uncontrolled diabetes.
- undergoing major dental work or an invasive medical procedure.

Safety and effectiveness have not been established for:

- patients under 18 years of age.
- pregnant or lactating women.
- patients with BMI ≥ 35 kg/m².
- single doses higher than 6 mL.

5. Summary of clinical evaluation and post-market clinical follow-up

5.1 Clinical background of the device

The gel on which Arthrosamid is based has been used for more than 20 years in products like Bulkamid[®] Hydrogel (used to treat urinary incontinence), and Arthramid[®] Vet (used for osteoarthritis in animals). Altogether, more than 1 million syringes of gels manufactured at Contura have been sold.

Five studies have been performed on Arthrosamid[®]. Additional evidence has been collected from scientific literature. The evidence shows that when Arthrosamid[®] is injected in the knee, it becomes part of the synovial membrane, a thin barrier that lines the inside of the joint (see figure on previous page). The treatment reduces pain and stiffness and improves physical function of the osteoarthritic knee, as well as quality of life. There is clinical evidence of a treatment effect up to 4 years after injection, but the expected total duration of effect is 5 years. The studies show that the treatment is safe.

5.2 The clinical evidence for CE-marking

The clinical studies on Arthrosamid® are described below.

A “proof of concept” study including 84 participants was conducted to obtain information on the effectiveness of Arthrosamid® in participants with osteoarthritis (OA) of the knee. In this study, a volume of 3 mL of Arthrosamid® was injected into the knee. In some participants, the injection was repeated once, after a minimum of 2 weeks (72 participants received a second injection). The study found that pain had significantly improved 4 months after the initial injection. Similarly, the stiffness and physical function were improved. The effect appeared to persist up to the last observation at 13 months. The study was published in 2018 [1].

Later, a retrospective study including 91 participants was conducted to obtain information on the safety of Arthrosamid® in patients with knee OA. Evaluation of the data showed no serious adverse events related to Arthrosamid® injection. No allergic reactions or infections were observed. Most of the adverse events reported by the participants can be described as normal symptoms of knee OA and were thus most probably not related to Arthrosamid®. An exception is the sensation of distention after treatment which was felt by 15 participants. Most (14 participants, corresponding to 93%) describe this sensation as short-lived (days to weeks). There were no indications of scarring after Arthrosamid® injection. This study was published in 2019 [2].

Another study was made on tissue removed from 7 patients, who had previously received treatment with Arthrosamid®, during total knee replacement. The overall finding of the investigation was that Arthrosamid® was found to be integrated in the synovial membrane, regardless of the time point after Arthrosamid® injection (5 – 33 months). No evidence of scar formation was seen.

The effectiveness and safety of one injection of 6 mL Arthrosamid® in patients with knee OA are evaluated in an ongoing study including 49 participants. Participants in this study have been followed for 4 years. Significant improvements have been seen in pain, stiffness and physical function at all evaluations, as well as a significant improvement in quality of life. A total of 14 adverse events related to Arthrosamid® have occurred in 12 participants since treatment. Of these, the most frequently reported events were joint pain, joint swelling and accumulation of fluid in the back of the knee. None of the events related to the use of Arthrosamid® were serious. Data evaluated at the 6- and 12-month time point were published in 2021 [4,5], 2-year data were published in 2022 [6], and 3-year data were published in 2023 [7,8]. Participants are invited to continue with yearly follow-up visits up to 5 years after treatment to evaluate safety and long-term performance of Arthrosamid®.

Another study is ongoing, which initially compared Arthrosamid® to Synvisc-One® for 1 year. Synvisc-One® is a different injectable (a hyaluronic acid) used for knee osteoarthritis. A total of 239 participants were randomised into the study. 120 participants received one injection (6 mL) of Synvisc-one® and 119 received one injection (6 mL) of Arthrosamid®. The study found no statistical difference in effectiveness between the two treatments. In the Arthrosamid® group, 35 participants experienced a total of 41 adverse events related to treatment. In the Synvisc-One® group, 9 participants experienced a total of 13 adverse events related to treatment. None of the related adverse events were serious. One-year data from this study were published in 2022 and 2023 [9-11]. After the first year, the study was extended to follow the participants in the Arthrosamid® group for up to 5 years. They have now been followed for 3 years, and significant improvements in pain, stiffness

and physical function have been found at all evaluations, as well as a significant improvement in quality of life.

5.3 Safety

Compared to total knee replacement, which is currently the go-to treatment when all other options have failed, treatment with Arthrosamid® is a minimally invasive, outpatient procedure. Arthrosamid® is neither absorbed nor degraded in the body, and data from 4 studies show a significant and long-lasting effect of the treatment. The collected evidence shows that most knee OA patients treated with Arthrosamid® had improvements in pain, stiffness and physical function of the OA knee, as well as quality of life. Most of the adverse events observed in the studies were short-lived and mild. The most frequent adverse events were joint pain and -swelling.

By 2023.DEC.31, a total of 5,098 patient packs (six-packs) were placed on the market. Safety data is continually collected from online databases, complaints, scientific literature and other sources. A review of safety data from launch to 2023.DEC.31 has found that Arthrosamid® is safe.

A post-market study evaluating the safety of Arthrosamid® was initiated in 2021 and is ongoing. In this study, 203 adult patients across 9 sites in Europe were treated with a single injection of 6 mL Arthrosamid®. After treatment with Arthrosamid®, participants will be followed for 5 years. No results are available yet.

6. Possible therapeutic alternatives

When considering alternative treatments, it is recommended that you contact your healthcare professional who will consider your individual situation.

Physicians may offer the following therapeutic alternatives to Arthrosamid®:

- Intra-articular treatments (treatments that are injected into the knee):
 - **Corticosteroids** may be effective in treating symptoms of inflammation, but the effect wears off quickly (2-6 weeks). Repeat injections are required for stable pain relief.
 - **Hyaluronic acid** has an effect lasting up to 6 months.
- Surgery: Surgeries that may be offered include total or partial joint replacement or knee osteotomy (a surgical procedure where a bone is cut and reshaped or realigned). Most patients have a good clinical outcome after surgery, but the risk of complications is relatively high compared to other treatments, and many patients will need revision surgery. However, surgery may be a good option for severely symptomatic patients in whom previous treatments have failed and who is experiencing a significant loss in quality of life.

The choice of treatment should be based on your specific symptoms, your goals and expectations, and the risks and benefits given your particular situation.

7. Suggested profile and training for users

Arthrosamid® is intended to be administered by qualified physicians familiar with intra-articular injection procedures, such as orthopaedic surgeons or rheumatologists.

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Původ	Č.	Rok	Účastník	Název
ISO	4180	2019	Obaly	Packaging- Complete, filled transport packages - General rules for the compilation of performance test schedules
DS/EN ISO	11607-1:2020 + A11:2022 + A1:2023	2020/2022/2023	Obaly	Packaging for terminally sterilised medical devices – Part 1: Requirements for materials, sterile barrier systems and packaging systems
EN ISO	11607-2:2020 + A11:2022 + A1:2023	2020/2022/2023	Obaly	Packaging for terminally sterilised medical devices – Part 2: Validation requirements for forming, sealing and assembly processes
ANSI/AAMI	ST67	2019	Sterilizace	Sterilization of health care products - requirements and guidance for selecting a sterility assurance level (SAL) for products labeled "sterile"
ASTM	F 1929-15	2015	Obaly	Standard Test Method for Detecting Seal Leaks in Porous Medical Packaging by Dye Penetration
ASTM	F 1980-21	2021	Obaly	Standard Guide for accelerated Aging of Sterile Barrier Systems and Medical Devices
ASTM	F 2096 – 11 (2019)	2019	Obaly	Standard Test Method for Detecting Gross Leaks in Medical Packaging by Internal Pressurisation (Bubble Test)
ASTM	F 1886/F1886M-16	2016	Obaly	Standard Test Method for Determining Integrity of Seals for Flexible Packaging by Visual Inspection
ASTM	F88/F 88M-23	2023	Obaly	Standard Test Method for Seal Strength of Flexible Barrier Materials
DS/EN ISO	13485:2016 + AC:2018 + A11:2021	2016/2018/2021	Hodnocení kvality – obecné	Medical devices – Quality management systems – Requirements for regulatory purposes
ISO	28594	2017	Postupy odběru vzorků	Combined accept-zero sampling systems and process control procedures for product acceptance
ISO	2859-1:1999 + Cor.1:2001 + Amd 1:2011	1999/2001/2011	Postupy odběru vzorků	Sampling procedures for inspection by attributes - Part 1: Sampling schemes indexed by acceptable quality limit (AQL) for lot-by-lot inspection
DS/EN/ ISO	11138-1	2017	Sterilizace	Sterilisation of health care products - biological indicators - Part 1
DS/EN ISO	11138-3	2017	Sterilizace	Sterilisation of health care products - biological indicators - Part 3 Biological indicators for moist heat sterilisation indicators
DS/EN ISO	11138-7	2019	Sterilizace	Sterilization of health care products - biological indicators – Part 7: Guidance for selection, use and interpretation of results
EN ISO	17665-1	2006	Sterilizace	Sterilisation of health care products – Moist heat – Part 1: Requirements for the development, validation and routine control of a sterilisation process for medical devices
EN	556-1:2001 + AC:2006	2001/2006	Sterilizace	Sterilisation of medical devices to be designated "STERILE" – Part 1: Requirements for terminally sterilised medical devices

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10. Literatura

1. Henriksen, M., Overgaard, A., Hartkopp, A., Bliddal, H., 2018. Intra-articular 2.5% polyacrylamide hydrogel for the treatment of knee osteoarthritis: an observational proof-of-concept cohort study. *Clinical and experimental rheumatology* 36, 1082-1085.
2. Overgaard, A., Bliddal, H., Henriksen, M., 2019. Safety of Intra-Articular Polyacrylamide Hydrogel for the Treatment of Knee Osteoarthritis Symptoms: A Retrospective Case Series. *Clin Ortho Adv Res J: COARJ-100001*.
3. Christensen, L., 2017. Case series: Report on a prospective histopathological study on biopsies from seven total arthroplasty knees (TKA) treated with polyacrylamide hydrogel (PAAG). Internal document.
4. Bliddal, H., Overgaard, A., Hartkopp, A., Beier, J., Conaghan, P.G., Henriksen, M., 2021. Polyacrylamide Hydrogel Injection for Knee Osteoarthritis: A 6 Months Prospective Study. *J Orthop Res Ther* 6, 1188.
5. Bliddal, H., Overgaard, A., Hartkopp, A., Beier, J., Conaghan, P.G., Henriksen, M., 2021. Polyacrylamide Hydrogel Injection for Knee Osteoarthritis: Results of a 52 Week Prospective Study. *Osteoarthritis and Cartilage* 29, S278.
6. Bliddal, H., Beier, J., Hartkopp, A., Conaghan, P.G., Henriksen, M., 2022. A Prospective Study of Polyacrylamide Hydrogel for Knee Osteoarthritis: Results From 2 Years After Treatment. *Osteoarthritis and Cartilage* 30, S371-S372.
7. Henriksen, M., Beier, J., Hartkopp, A., Conaghan, P.G., Bliddal, H., 2023. 3 year results from a prospective study of polyacrylamide hydrogel for knee osteoarthritis. *Osteoarthritis and Cartilage* 31, P682-3.
8. Bliddal, H., Beier, J., Hartkopp, A., Conaghan, P.G., Henriksen, M., 2023. AB1185 Polyacrylamide hydrogel for the treatment of knee osteoarthritis: 3 year follow up results of a prospective clinical study. *Annals of the Rheumatic Diseases* 82:1822-1823.
9. Bliddal, H., Beier, J., Hartkopp, A., Conaghan, P.G., Henriksen, M., 2022. One-Year Performance of Polyacrylamide Hydrogel vs. Hyaluronic Acid: A Randomised Controlled Study. *Osteoarthritis and Cartilage* 30, S370-S371.
10. Bliddal, H., Beier, J., Hartkopp, A., Conaghan, P.G., Henriksen, M., 2022. One-Year Performance of Polyacrylamide Hydrogel vs. Hyaluronic Acid in Age, BMI and Kellgren-

Lawrence Subgroups: a Subgroup Analysis of a Randomised Study. *Osteoarthritis and Cartilage* 30, S373-S374.

11. Bliddal, H., Beier, J., Hartkopp, A., Conaghan, P.G., Henriksen, M., 2023. Polyacrylamide hydrogel versus hyaluronic acid in knee osteoarthritis: a subgroup analysis of a randomized study. *Orthopaedic Proceedings* 105-B(SUPP 13):81.
12. Christensen, L., Daugaard, S., 2016. Histological Appearance of the Synovial Membrane after Treatment of Knee Osteoarthritis with Polyacrylamide Gel Injections: A Case Report. *Journal of Arthritis* 5, 217.
13. Maulana, R., Cole, A., Lee, P.Y.F., 2022. Reduction in Patellofemoral Bone Marrow Lesions Following Single Arthrosamid Intra-Articular Injection of Polyacrylamide Hydrogel (ipaag) in the Treatment of Advanced Osteoarthritis. *J Arthritis* 11, 24-26.
14. Negredo, E., Puig, J., Ornelas, A., Echeverria, P., Bonjoch, A., Estany, C., Higuera, C., Gonzalez-Mestre, V., Clotet, B., 2015. Ten-Year Safety with Polyacrylamide Gel Used to Correct Facial Lipoatrophy in HIV-Infected Patients. *AIDS Res Hum Retroviruses* 31, 817-821.
15. Mouritsen, L., Lose, G., Moller-Bek, K., 2014. Long-term follow-up after urethral injection with polyacrylamide hydrogel for female stress incontinence. *Acta Obstet Gynecol Scand* 93, 209-212.
16. Bruyere, O., Honvo, G., Veronese, N., Arden, N.K., Branco, J., Curtis, E.M., Al-Daghri, N.M., Herrero-Beaumont, G., Martel-Pelletier, J., Pelletier, J.P., Rannou, F., Rizzoli, R., Roth, R., Uebelhart, D., Cooper, C., Reginster, J.Y., 2019. An updated algorithm recommendation for the management of knee osteoarthritis from the European Society for Clinical and Economic Aspects of Osteoporosis, Osteoarthritis and Musculoskeletal Diseases (ESCEO). *Seminars in arthritis and rheumatism* 49, 337-350.
17. Najm, A., Alunno, A., Gwinnutt, J.M., Weill, C., Berenbaum, F., 2021. Efficacy of intra-articular corticosteroid injections in knee osteoarthritis: A systematic review and meta-analysis of randomized controlled trials. *Joint Bone Spine* 88, 105198.
18. Tammachote, N., Kanitnate, S., Yakumpor, T., Panichkul, P., 2016. Intra-Articular, Single-Shot Hylan G-F 20 Hyaluronic Acid Injection Compared with Corticosteroid in Knee Osteoarthritis: A Double-Blind, Randomized Controlled Trial. *The Journal of bone and joint surgery. American volume* 98, 885-892.
19. Wang, C.P., Lee, W.C., Hsieh, R.L., 2022. Effects of Repeated Co-Injections of Corticosteroids and Hyaluronic Acid on Knee Osteoarthritis: A Prospective, Double-Blind Randomized Controlled Trial. *Am J Med.* 135, 641-649.
20. Marsland, D., Mumith, A., Barlow, I.W., 2014. Systematic review: the safety of intra-articular corticosteroid injection prior to total knee arthroplasty. *The Knee* 21, 6-11.

21. Gomis, A., Pawlak, M., Balazs, E.A., Schmidt, R.F., Belmonte, C., 2004. Effects of different molecular weight elastoviscous hyaluronan solutions on articular nociceptive afferents. *Arthritis and rheumatism* 50, 314-326.
22. Altman, R.D., Bedi, A., Karlsson, J., Sancheti, P., Schemitsch, E., 2016. Product Differences in Intra-articular Hyaluronic Acids for Osteoarthritis of the Knee. *The American journal of sports medicine* 44, 2158-2165.
23. Sun, S.F., Hsu, C.W., Lin, H.S., Liou, I.H., Chen, Y.H., Hung, C.L., 2017. Comparison of Single Intra-Articular Injection of Novel Hyaluronan (HYA-JOINT Plus) with Synvisc-One for Knee Osteoarthritis: A Randomized, Controlled, Double-Blind Trial of Efficacy and Safety. *The Journal of bone and joint surgery. American volume* 99, 462-471.
24. Zhao, H., Liu, H., Liang, X., Li, Y., Wang, J., Liu, C., 2016. Hylan G-F 20 Versus Low Molecular Weight Hyaluronic Acids for Knee Osteoarthritis: A Meta-Analysis. *BioDrugs : clinical immunotherapeutics, biopharmaceuticals and gene therapy* 30, 387-396.
25. Balazs, E.A., 2004. Viscosupplementation for treatment of osteoarthritis: from initial discovery to current status and results. *Surgical technology international* 12, 278-289.
26. He, W.W., Kuang, M.J., Zhao, J., Sun, L., Lu, B., Wang, Y., Ma, J.X., Ma, X.L., 2017. Efficacy and safety of intraarticular hyaluronic acid and corticosteroid for knee osteoarthritis: A meta-analysis. *International journal of surgery (London, England)* 39, 95-103.
27. Leighton, R., Akemark, C., Therrien, R., Richardson, J.B., Andersson, M., Todman, M.G., Arden, N.K., 2014. NASHA hyaluronic acid vs. methylprednisolone for knee osteoarthritis: a prospective, multi-centre, randomized, non-inferiority trial. *Osteoarthritis and cartilage / OARS, Osteoarthritis Research Society* 22, 17-25.
28. Phillips, M., Vannabouathong, C., Devji, T., Patel, R., Gomes, Z., Patel, A., Dixon, M., Bhandari, M., 2020. Differentiating factors of intra-articular injectables have a meaningful impact on knee osteoarthritis outcomes: a network meta-analysis. *Knee surgery, sports traumatology, arthroscopy: official journal of the ESSKA* 28, 3031-3039.
29. Leighton, R., Fitzpatrick, J., Smith, H., Crandall, D., Flannery, C.R., Conrozier, T., 2018. Systematic clinical evidence review of NASHA (Durolane hyaluronic acid) for the treatment of knee osteoarthritis. *Open Access Rheumatol* 10, 43-54.
30. Maheu, E., Bannuru, R.R., Herrero-Beaumont, G., Allali, F., Bard, H., Migliore, A., 2019. Why we should definitely include intra-articular hyaluronic acid as a therapeutic option in the management of knee osteoarthritis: Results of an extensive critical literature review. *Seminars in arthritis and rheumatism* 48, 563-572.
31. Maheu, E., Rannou, F., Reginster, J.Y., 2016. Efficacy and safety of hyaluronic acid in the management of osteoarthritis: Evidence from real-life setting trials and surveys. *Seminars in arthritis and rheumatism* 45, S28-33.

32. Tiendrebego, E., Choueiri, M., Chevalier, X., Conrozier, T., Eymard, F., 2021. Does the Presence of Neuropathic Pain Influence the Response to Hyaluronic Acid in Patients with Knee Osteoarthritis? *Cartilage* 13, 1548S-1556S.
33. Singh, H., Knapik, D.M., Polce, E.M., Eikani, C.K., Bjornstad, A.H., Gursoy, S., Perry, A.K., Westrick, J.C., Yanke, A.B., Verma, N.N., Cole, B.J., Chahla, J.A., 2022. Relative Efficacy of Intra-articular Injections in the Treatment of Knee Osteoarthritis: A Systematic Review and Network Meta-analysis. *Am J Sports Med* 50, 3140-3148.
34. Jevsevar, D., Donnelly, P., Brown, G.A., Cummins, D.S., 2015. Viscosupplementation for Osteoarthritis of the Knee: A Systematic Review of the Evidence. *The Journal of bone and joint surgery. American volume* 97, 2047-2060.
35. Ke, Y., Jiang, W., Xu, Y., Chen, Y., Zhang, Q., Xue, Q., Lin, J., Ngai, W., Nian, G., Fazeli, M.S., Xie, Y., Zhu, Z., 2021. Efficacy and safety of a single intra-articular injection of 6 ml Hylan G-F 20 compared to placebo in Chinese patients with symptomatic knee osteoarthritis : C-SOUND study, a 26-week multicenter double-blind randomized placebo-controlled trial in China. *BMC Musculoskelet Disord* 22, 428.
36. Pereira, T.V., Jüni, P., Saadat, P., Xing, D., Yao, L., Bobos, P., Agarwal, A., Hincapié, C.A., da Costa, B.R., 2022. Viscosupplementation for knee osteoarthritis: systematic review and meta-analysis. *BMJ* 378, e069722.
37. Cui, B., Chen, Y., Tian, Y., Liu, H., Huang, Y., Yin, G., Xie, Q., 2022. Effects of medications on incidence and risk of knee and hip joint replacement in patients with osteoarthritis: a systematic review and meta-analysis. *Adv Rheumatol* 62, 22.
38. Miller, L.E., Block, J.E., 2013. US-Approved Intra-Articular Hyaluronic Acid Injections are Safe and Effective in Patients with Knee Osteoarthritis: Systematic Review and Meta-Analysis of Randomized, Saline-Controlled Trials. *Clinical medicine insights. Arthritis and musculoskeletal disorders* 6, 57-63.
39. Petterson, S.C., Plancher, K.D., 2019. Single intra-articular injection of lightly cross-linked hyaluronic acid reduces knee pain in symptomatic knee osteoarthritis: a multicenter, double-blind, randomized, placebo-controlled trial. *Knee surgery, sports traumatology, arthroscopy : official journal of the ESSKA* 27, 1992-2002.
40. Miller, L.E., Bhattacharyya, S., Parrish, W.R., Fredericson, M., Bisson, B., Altman, R.D., 2021. Safety of Intra-Articular Hyaluronic Acid for Knee Osteoarthritis: Systematic Review and Meta-Analysis of Randomized Trials Involving More than 8,000 Patients. *Cartilage* 13, 351S-363S.
41. Bagga, H., Burkhardt, D., Sambrook, P., March, L., 2006. Longterm effects of intraarticular hyaluronan on synovial fluid in osteoarthritis of the knee. *J Rheumatol* 33, 946-950.
42. Xiao, J., Hu, Y., Huang, L., Huang, Z.F., Jiang, W.Z., Luo, Y.Q., Jia, M.Y., Chen, D., Shi, Z.J., 2021. Injection route affects intra-articular hyaluronic acid distribution and clinical outcome in viscosupplementation treatment for knee osteoarthritis: a combined cadaver study and randomized clinical trial. *Drug Deliv Transl Res* 11, 279-291.

43. Benazzo, F., Perticarini, L., Padolino, A., Castelli, A., Gifuni, P., Lovato, M., Manzini, C., Giordan, N., 2016. A multi-centre, open label, long-term follow-up study to evaluate the benefits of a new viscoelastic hydrogel (Hymovis(R)) in the treatment of knee osteoarthritis. *European review for medical and pharmacological sciences* 20, 959-968.
44. Bisicchia, S., Bernardi, G., Tudisco, C., 2016. HYADD 4 versus methylprednisolone acetate in symptomatic knee osteoarthritis: a single-centre single blind prospective randomised controlled clinical study with 1-year follow-up. *Clinical and experimental rheumatology* 34, 857-863.
45. Hermans, J., Bierma-Zeinstra, S.M.A., Bos, P.K., Niesten, D.D., Verhaar, J.A.N., Reijman, M., 2019. The effectiveness of high molecular weight hyaluronic acid for knee osteoarthritis in patients in the working age: a randomised controlled trial. *BMC musculoskeletal disorders* 20, 196.
46. Yan, C.H., Chan, W.L., Yuen, W.H., Yung, P.S., Ip, K.Y., Fan, J.C., Chiu, K.Y., 2015. Efficacy and safety of hylan G-F 20 injection in treatment of knee osteoarthritis in Chinese patients: results of a prospective, multicentre, longitudinal study. *Hong Kong medical journal = Xianggang yi xue za zhi* 21, 327-332.
47. Price, A.J., Alvand, A., Troelsen, A., Katz, J.N., Hooper, G., Gray, A., Carr, A., Beard, D., 2018. Knee replacement. *Lancet* 392, 1672-1682.
48. Peng G, Liu M, Guan Z, Hou Y, Liu Q, Sun X, Zhu X, Feng W, Zeng J, Zhong Z, Zeng Y. Patellofemoral arthroplasty versus total knee arthroplasty for isolated patellofemoral osteoarthritis: a systematic review and meta-analysis. *J Orthop Surg Res.* 2021 Apr 15;16(1):264.
49. Liddle, A.D., Judge, A., Pandit, H., Murray, D.W., 2014. Adverse outcomes after total and unicompartmental knee replacement in 101,330 matched patients: a study of data from the National Joint Registry for England and Wales. *Lancet* 384, 1437-1445.
50. Li C, Li Z, Shi L, Gao F, Sun W. The short-term effectiveness and safety of second-generation patellofemoral arthroplasty and total knee arthroplasty on isolated patellofemoral osteoarthritis: a systematic review and meta-analysis. *J Orthop Surg Res.* 2021 Jun 2;16(1):358.
51. Olsen, U., Lindberg, M.F., Rose, C., Denison, E., Gay, C., Aamodt, A., Brox, J.I., Skare, Ø., Furnes, O., Lee, K., Lerdal, A., 2022. Factors Correlated With Physical Function 1 Year After Total Knee Arthroplasty in Patients With Knee Osteoarthritis: A Systematic Review and Meta-analysis. *JAMA Netw Open* 5, e2219636.
52. Carr, A.J., Robertsson, O., Graves, S., Price, A.J., Arden, N.K., Judge, A., Beard, D.J., 2012. Knee replacement. *Lancet* 379, 1331-1340.
53. Belmont, P.J., Jr., Goodman, G.P., Waterman, B.R., Bader, J.O., Schoenfeld, A.J., 2014. Thirty-day postoperative complications and mortality following total knee arthroplasty: incidence and risk factors among a national sample of 15,321 patients. *The Journal of bone and joint surgery. American volume* 96, 20-26.

54. Wilson, M.G., Kelley, K., Thornhill, T.S., 1990. Infection as a complication of total knee-replacement arthroplasty. Risk factors and treatment in sixty-seven cases. The Journal of bone and joint surgery. American volume 72, 878-883.
55. Skou, S.T., Roos, E.M., Laursen, M.B., Rathleff, M.S., Arendt-Nielsen, L., Simonsen, O., Rasmussen, S., 2015. A Randomized, Controlled Trial of Total Knee Replacement. The New England journal of medicine 373, 1597-1606.
56. Skou, S.T., Roos, E.M., Laursen, M.B., Rathleff, M.S., Arendt-Nielsen, L., Rasmussen, S., Simonsen, O., 2018. Total knee replacement and non-surgical treatment of knee osteoarthritis: 2-year outcome from two parallel randomized controlled trials. Osteoarthritis and cartilage / OARS, Osteoarthritis Research Society 26, 1170-1180.

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Quelle	Nr.	Jahr	Thema	Titel
DS/EN ISO	10993-1	2020	Biologische Untersuchung	Biological evaluation of medical devices – Part 1 Evaluation and testing
DS/EN ISO	10993-3	2014	Biologische Untersuchung	Biological evaluation of medical devices – Part 3: Test for genotoxicity, carcinogenicity and reproductive toxicity
DS/EN ISO	10993-5	2009	Biologische Untersuchung	Biological evaluation of medical devices – Part 5: Tests for in vitro cytotoxicity
DS/EN ISO	10993-6	2016	Biologische Untersuchung	Biological evaluation of medical devices – Part 6: Tests for local effects after implantation
DS/EN ISO	10993-9	2021	Biologische Untersuchung	Biological evaluation of medical devices – Part 9 Framework for identification and degradation of potential degradation products
DS/EN ISO	10993-10	2023	Biologische Untersuchung	Biological evaluation of medical devices – Part 10: Tests for irritation and delayed-type hypersensitivity
DS/EN ISO	10993-11	2018	Biologische Untersuchung	Biological evaluation of medical devices – Part 11: Tests for systemic toxicity
DS/EN ISO	10993-12	2021	Biologische Untersuchung	Biological evaluation of medical devices - Part 12: Sample preparation and reference materials
DS/EN ISO	10993-13	2010	Biologische Untersuchung	Biological evaluation of medical devices – Part 13: Identification and quantification of degradation products from polymeric medicinal devices
DS/EN ISO	10993-18:2020 + A1:2023	2020/2023	Biologische Untersuchung	Biological evaluation of medical devices – Part 18: Chemical characterisation of materials
DS/ISO/TS	10993-19	2020	Biologische Untersuchung	Biological evaluation of medical devices – Part 19: Physico-chemical, morphological and topographical characterisation of materials
DS/EN/TR	10993-33	2015	Biologische Untersuchung	Biological evaluation of medical devices - Part 33: Guidance on tests to evaluate genotoxicity - Supplement to ISO 10993-3
DS/EN/TR	21726	2019	Biologische Untersuchung	Biological evaluation of medical devices - Application of the threshold of toxicological concern (TTC) for assessing biocompatibility of medical devices constituents
DS/EN	12469	2000	Reinräume	Biotechnology - Performance criteria for microbiological safety cabinets
EN ISO	14644-1	2015	Reinräume	Cleanrooms and associated controlled environments – Part 1: Classification of air cleanliness
EN ISO	14644-2	2015	Reinräume	Cleanrooms and Associated Controlled Environments - Part 2: Specification for Testing and Monitoring to Prove Continued Compliance with ISO 14644-1
EN ISO	14644-3:2019 + Cor. 2020-06	2019/2020	Reinräume	Cleanrooms and Associated Controlled Environments - Part 3: Test methods
EN ISO	14644-4	2022	Reinräume	Cleanrooms and Associated Controlled Environments - Part 4: Design, Construction and Start-up
EN ISO	14644-5	2004	Reinräume	Cleanrooms and Associated Controlled Environments - Part 5: Entwurf
EN ISO	14644-7	2004	Reinräume	Cleanrooms and Associated Controlled Environments - Part 7: Separative devices (clean air hood, gloveboxes, isolators and mini-environments)
EN ISO	14644-8	2022	Reinräume	Cleanrooms and Associated Controlled Environments - Part 8: Classification of airborne molecular contamination
DS/EN	17141	2020	Reinräume	Cleanrooms and associated controlled environments - Bio-contamination control

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Quelle	Nr.	Jahr	Thema	Titel
Ph. Eur.	2.2.44	2024	Herstellung und Validierung	Total Organic Carbon (Gesamtgehalt an organischem Kohlenstoff) (TOC)
Ph. Eur.	2.2.6	2024	Herstellung und Validierung	Refractive Index (Refraktionsindex)
Ph. Eur.	2.4.8	2024	Herstellung und Validierung	Heavy metals (Schwermetalle)
Ph. Eur.	2.6.12	2024	Herstellung und Validierung	Microbiological Examination of Non-Sterile Products: microbial enumeration tests
Ph. Eur.	1238	2024	Herstellung und Validierung	Air Medicinal (Medizinische Luft)
DS/ISO	2859-1:1999 + Amd.1:2011	1999/2011	Herstellung und Validierung	Sampling procedures for inspection by attributes • Part 1: Sampling schemes indexed by acceptable quality limit (AQL) for lot-by-lot inspection
ISO	4180	2019	Verpackung	Packaging- Complete, filled transport packages - General rules for the compilation of performance test schedules
DS/EN ISO	11607-1:2020 + A11:2022 + A1:2023	2020/2022/2023	Verpackung	Packaging for terminally sterilised medical devices – Part 1: Requirements for materials, sterile barrier systems and packaging systems
EN ISO	11607-2:2020 + A11:2022 + A1:2023	2020/2022/2023	Verpackung	Packaging for terminally sterilised medical devices – Part 2: Validation requirements for forming, sealing and assembly processes
ANSI/AAMI	ST67	2019	Sterilisation	Sterilization of health care products - requirements and guidance for selecting a sterility assurance level (SAL) for products labeled "sterile"
ASTM	F 1929-15	2015	Verpackung	Standard Test Method for Detecting Seal Leaks in Porous Medical Packaging by Dye Penetration
ASTM	F 1980-21	2021	Verpackung	Standard Guide for accelerated Aging of Sterile Barrier Systems and Medical Devices
ASTM	F 2096 – 11 (2019)	2019	Verpackung	Standard Test Method for Detecting Gross Leaks in Medical Packaging by Internal Pressurisation (Bubble Test)
ASTM	F 1886/F1886M-16	2016	Verpackung	Standard Test Method for Determining Integrity of Seals for Flexible Packaging by Visual Inspection
ASTM	F88 88M-23	2023	Verpackung	Standard Test Method for Seal Strength of Flexible Barrier Materials
DS/EN ISO	13485:2016 + AC:2018 + A11:2021	2016/2018/2021	QS-Allgemein	Medical devices – Quality management systems – Requirements for regulatory purposes
ISO	28594	2017	Stichprobenverfahren	Combined accept-zero sampling systems and process control procedures for product acceptance
ISO	2859-1:1999 + Cor.1:2001 + Amd 1:2011	1999/2001/2011	Stichprobenverfahren	Sampling procedures for inspection by attributes - Part 1: Sampling schemes indexed by acceptable quality limit (AQL) for lot-by-lot inspection
DS/EN/ ISO	11138-1	2017	Sterilisation	Sterilisation of health care products - biological indicators - Part 1
DS/EN ISO	11138-3	2017	Sterilisation	Sterilisation of health care products - biological indicators - Part 3 Biological indicators for moist heat sterilisation indicators
DS/EN ISO	11138-7	2019	Sterilisation	Sterilization of health care products - biological indicators – Part 7: Guidance for selection, use and interpretation of results
EN ISO	17665-1	2006	Sterilisation	Sterilisation of health care products – Moist heat – Part 1: Requirements for the development, validation and routine control of a sterilisation process for medical devices
EN	556-1:2001 + AC:2006	2001/2006	Sterilisation	Sterilisation of medical devices to be designated “STERILE” – Part 1: Requirements for terminally sterilised medical devices

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Alkuperä	Nro	Vuosi	Aihe	Nimi
Ph.Eur	2.2.3	2024	Valmistus ja validointi	Potentiometric Determination of pH
Ph.Eur	2.2.32	2024	Valmistus ja validointi	Loss on Drying
Ph.Eur	2.2.44	2024	Valmistus ja validointi	Total Organic Carbon (TOC)
Ph.Eur	2.2.6	2024	Valmistus ja validointi	Refractive Index
Ph.Eur	2.4.8	2024	Valmistus ja validointi	Heavy metals
Ph.Eur	2.6.12	2024	Valmistus ja validointi	Microbiological Examination of Non-Sterile Products: microbial enumeration tests
Ph.eur	1238	2024	Valmistus ja validointi	Air Medicinal
DS/ISO	2859-1:1999 + Amd.1:2011	1999/2011	Valmistus ja validointi	Sampling procedures for inspection by attributes • Part 1: Sampling schemes indexed by acceptable quality limit (AQL) for lot-by-lot inspection
ISO	4180	2019	Pakkaaminen	Packaging- Complete, filled transport packages - General rules for the compilation of performance test schedules
DS/EN ISO	11607-1:2020 + A11:2022 + A1:2023	2020/2022/2023	Pakkaaminen	Packaging for terminally sterilised medical devices – Part 1: Requirements for materials, sterile barrier systems and packaging systems
EN ISO	11607-2:2020 + A11:2022 + A1:2023	2020/2022/2023	Pakkaaminen	Packaging for terminally sterilised medical devices – Part 2: Validation requirements for forming, sealing and assembly processes
ANSI/AAMI	ST67	2019	Sterilointi	Sterilization of health care products - requirements and guidance for selecting a sterility assurance level (SAL) for products labeled "sterile"
ASTM	F 1929-15	2015	Pakkaaminen	Standard Test Method for Detecting Seal Leaks in Porous Medical Packaging by Dye Penetration
ASTM	F 1980-21	2021	Pakkaaminen	Standard Guide for accelerated Aging of Sterile Barrier Systems and Medical Devices
ASTM	F 2096 – 11 (2019)	2019	Pakkaaminen	Standard Test Method for Detecting Gross Leaks in Medical Packaging by Internal Pressurisation (Bubble Test)
ASTM	F 1886/F1886M-16	2016	Pakkaaminen	Standard Test Method for Determining Integrity of Seals for Flexible Packaging by Visual Inspection
ASTM	F88/F 88M-23	2023	Pakkaaminen	Standard Test Method for Seal Strength of Flexible Barrier Materials
DS/EN ISO	13485:2016 + AC:2018 + A11:2021	2016/2018/2021	Laadunarviointi – yleinen	Terveysthuollon laitteet ja tarvikkeet. Laadunhallintajärjestelmät. Vaatimukset viranomaismääräyksiä varten.
ISO	28594	2017	Näytteenottomenetelmät	Combined accept-zero sampling systems and process control procedures for product acceptance
ISO	2859-1:1999 + Korj.1.:2001 + Lis. 1:2011	1999/2001/2011	Näytteenottomenetelmät	Sampling procedures for inspection by attributes - Part 1: Sampling schemes indexed by acceptable quality limit (AQL) for lot-by-lot inspection
DS/EN/ ISO	11138-1	2017	Sterilointi	Sterilisation of health care products - biological indicators - Part 1
DS/EN ISO	11138-3	2017	Sterilointi	Sterilisation of health care products - biological indicators - Part 3 Biological indicators for moist heat sterilisation indicators
DS/EN ISO	11138-7	2019	Sterilointi	Sterilization of health care products - biological indicators – Part 7: Guidance for selection, use and interpretation of results
EN ISO	17665-1	2006	Sterilointi	Sterilisation of health care products – Moist heat – Part 1: Requirements for the development, validation and routine control of a sterilisation process for medical devices
EN	556-1:2001 + AC:2006	2001/2006	Sterilointi	Sterilisation of medical devices to be designated "STERILE" – Part 1: Requirements for terminally sterilised medical devices

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5.2.2	Safety assessment study	15
5.2.2.1	Clinical investigation design	15
5.2.2.2	Results	17
5.2.2.3	Conclusion	18
5.2.3	Histopathological report	19
5.2.3.1	Clinical investigation design	19
5.2.3.2	Results	19
5.2.3.3	Conclusion	22
5.2.4	Prospective, open-label, clinical investigation of intra-articular polyacrylamide hydrogel injection in subjects with knee OA ("IDA") (CIV ID: CIV-19-07-029146)	22
5.2.4.1	Clinical investigation design	23
5.2.4.2	Results	25
5.2.4.3	Conclusion	34
5.2.5	A multi-centre randomised, controlled, double-blind clinical investigation of intra-articular polyacrylamide hydrogel injection in subjects with knee osteoarthritis ("ROSA") (CIV ID: CIV-19-01-026964)	34
5.2.5.1	Clinical investigation design	35
5.2.5.2	Results	37
5.2.5.3	Conclusion	46
5.3	Other clinical data	47
5.4	Summary of the clinical performance and safety	47
5.5	Post-market clinical follow-up	48
6.	Possible therapeutic alternatives	49
6.1	Intra-articular treatments	49
6.1.1	Corticosteroids	49
6.1.2	Hyaluronic acid	49
6.1.3	Noltrex	50
6.2	Surgery	50
7.	Suggested profile and training for users	51
8.	Harmonised standards and applied Common Specifications	51
9.	Revision history	54
10.	Bibliography	55
	Samantekt á öryggi og klínískri virkni - hluti II – fyrir sjúklinga	61
1.	Auðkenning tækis og almennar upplýsingar	61
1.1	Vöruheiti tækis	61
1.2	Framleiðandi	61

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Study endpoints		Study results
Primary endpoint	To evaluate the change from baseline in pain (WOMAC) at 4 months after injection of “Hydrogel B” into the knee joint.	A statistically significant improvement in pain at 4 months after injection of “Hydrogel B” into the knee joint was demonstrated (p<0.0001).
Secondary endpoints	To evaluate changes from baseline in the physical function and stiffness sub-scales of WOMAC as well as the total WOMAC score at 4, 7, and 13 months.	Improvement was seen in the other WOMAC domains (stiffness, function, and total) after 4 months. The effect seems to be long-term with effectiveness maintained for 13 months.

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The analysis of the change from baseline in the transformed WOMAC pain subscale (0-100) showed a significant reduction ($p < 0.0001$) of 22.0 points from baseline to year 4 for the full ITT analysis set and 21.8 points for the year-2 extension subset, indicating reduced pain in the treated knee. The most conservative approach, BOCF, also resulted in a significant reduction ($p < 0.0001$) of 13.0 points from baseline to year 4.

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WOMAC stiffness subscale

The scores for the transformed WOMAC stiffness subscale at baseline, week 13, week 26, week 52, year 2, year 3 and year 4 are summarised for the ITT analysis set together with the respective changes from baseline in [Table 11](#). The mean changes from baseline are shown graphically in [Mynd](#)

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week 26, week 52, year 2, year 3 and year 4 are summarised for the ITT analysis set together with the respective changes from baseline in [Table 13](#). The changes from baseline are shown graphically in [Figure 6](#) for the ITT analysis set, both as observed and using the BOCF approach, and for the year-2 extension subset.

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WOMAC stiffness subscale

The analysis of the change from baseline in the transformed WOMAC stiffness subscale (0-100) was performed using the same MMRM as for the primary endpoint. The analysis showed a statistically significant treatment difference at week 4 in favour of Synvisc-one ($p = 0.0315$) while treatment differences in favour of Arthrosamid® at weeks 12, 26 and 52 were not significant ([Table 17](#)).

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OMERACT-OARSI response

The analysis of the proportion of positive responders to the OMERACT-OARSI criteria was performed using a repeated measures logistic model treatment. The analysis showed a treatment ratio (odds ratio) in favour of Synvisc-one at week 4 that was not statistically significant ([Table 20](#)). At weeks 12, 26 and 52, the treatment ratio was close to 1, indicating no difference between the 2 groups.

Table 20. Analysis of OMERACT-OARSI criteria - ITT analysis set.

	E/n	(%)	Synvisc-One / PAAG-OA Odds ratio (95% CI)	p-value
Week 4				
Synvisc-One	74 / 120	(61.7)		
PAAG-OA	64 / 119	(53.8)	1.38 (0.82; 2.31)	0.2200
Week 12				
Synvisc-One	77 / 119	(64.7)		
PAAG-OA	79 / 115	(68.7)	0.87 (0.50; 1.49)	0.5994
Week 26				
Synvisc-One	76 / 117	(65.0)		
PAAG-OA	76 / 114	(66.7)	0.95 (0.55; 1.62)	0.8466
Week 52				
Synvisc-One	62 / 109	(56.9)		
PAAG-OA	60 / 107	(56.1)	1.01 (0.59; 1.73)	0.9612

CI: confidence interval; E: Number of OMERACT-OARSI responders; n: number with non-missing values.

EuroQoL-5D-5L

The answers from each of the 5 dimensions from EQ-5D-5L were combined into one measurement termed the EQ-5D-5L. The analysis of the change from baseline in the EQ-5D-5L health state showed a statistically significant treatment difference in favour of Synvisc-one at week 4 while the results at weeks 12, 26 and 52 were similar between the Arthroamid® arm and the Synvisc-one arm ([Table 21](#)).

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6.1.2 Hyaluronic acid

There are numerous HA-based viscosupplemental products for human use on the market. The preparations have varying average molecular weight (0.5 6.0 million Da) and therefore show great variations in visco-elastic properties, in turn leading to varying effects in terms of pain reduction and functional improvement [21-24] as well as required dosage and frequency of injections [25].

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12. Christensen, L., Daugaard, S., 2016. Histological Appearance of the Synovial Membrane after Treatment of Knee Osteoarthritis with Polyacrylamide Gel Injections: A Case Report. *Journal of Arthritis* 5, 217.
13. Maulana, R., Cole, A., Lee, P.Y.F., 2022. Reduction in Patellofemoral Bone Marrow Lesions Following Single Arthrosamid Intra -Articular Injection of Polyacrylamide Hydrogel (ipaag) in the Treatment of Advanced Osteoarthritis. *J Arthritis* 11, 24-26.
14. Negredo, E., Puig, J., Ornelas, A., Echeverria, P., Bonjoch, A., Estany, C., Higuera, C., Gonzalez-Mestre, V., Clotet, B., 2015. Ten-Year Safety with Polyacrylamide Gel Used to Correct Facial Lipoatrophy in HIV-Infected Patients. *AIDS Res Hum Retroviruses* 31, 817-821.
15. Mouritsen, L., Lose, G., Moller-Bek, K., 2014. Long-term follow-up after urethral injection with polyacrylamide hydrogel for female stress incontinence. *Acta Obstet Gynecol Scand* 93, 209-212.
16. Bruyere, O., Honvo, G., Veronese, N., Arden, N.K., Branco, J., Curtis, E.M., Al-Daghri, N.M., Herrero-Beaumont, G., Martel-Pelletier, J., Pelletier, J.P., Rannou, F., Rizzoli, R., Roth, R., Uebelhart, D., Cooper, C., Reginster, J.Y., 2019. An updated algorithm recommendation for the management of knee osteoarthritis from the European Society for Clinical and Economic Aspects of Osteoporosis, Osteoarthritis and Musculoskeletal Diseases (ESCEO). *Seminars in arthritis and rheumatism* 49, 337-350.
17. Najm, A., Alunno, A., Gwinnutt, J.M., Weill, C., Berenbaum, F., 2021. Efficacy of intra-articular corticosteroid injections in knee osteoarthritis: A systematic review and meta-analysis of randomized controlled trials. *Joint Bone Spine* 88, 105198.
18. Tammachote, N., Kanitnate, S., Yakumpor, T., Panichkul, P., 2016. Intra-Articular, Single-Shot Hylan G-F 20 Hyaluronic Acid Injection Compared with Corticosteroid in Knee Osteoarthritis: A Double-Blind, Randomized Controlled Trial. *The Journal of bone and joint surgery. American volume* 98, 885-892.
19. Wang, C.P., Lee, W.C., Hsieh, R.L., 2022. Effects of Repeated Co-Injections of Corticosteroids and Hyaluronic Acid on Knee Osteoarthritis: A Prospective, Double-Blind Randomized Controlled Trial. *Am J Med.* 135, 641-649.
20. Marsland, D., Mumith, A., Barlow, I.W., 2014. Systematic review: the safety of intra-articular corticosteroid injection prior to total knee arthroplasty. *The Knee* 21, 6-11.
21. Gomis, A., Pawlak, M., Balazs, E.A., Schmidt, R.F., Belmonte, C., 2004. Effects of different molecular weight elastoviscous hyaluronan solutions on articular nociceptive afferents. *Arthritis and rheumatism* 50, 314-326.
22. Altman, R.D., Bedi, A., Karlsson, J., Sancheti, P., Schemitsch, E., 2016. Product Differences in Intra-articular Hyaluronic Acids for Osteoarthritis of the Knee. *The American journal of sports medicine* 44, 2158-2165.

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Origine	N.	Anno	Oggetto	Titolo
DS/EN ISO	10993-1	2020	Valutazione biologica	Biological evaluation of medical devices – Part 1: Evaluation and testing
DS/EN ISO	10993-3	2014	Valutazione biologica	Biological evaluation of medical devices – Part 3: Test for genotoxicity, carcinogenicity and reproductive toxicity

Origine	N.	Anno	Oggetto	Titolo
DS/EN ISO	10993-5	2009	Valutazione biologica	Biological evaluation of medical devices – Part 5: Tests for in vitro cytotoxicity
DS/EN ISO	10993-6	2016	Valutazione biologica	Biological evaluation of medical devices – Part 6: Tests for local effects after implantation
DS/EN ISO	10993-9	2021	Valutazione biologica	Biological evaluation of medical devices – Part 9: Framework for identification and degradation of potential degradation products
DS/EN ISO	10993-10	2023	Valutazione biologica	Biological evaluation of medical devices – Part 10: Tests for irritation and delayed-type hypersensitivity
DS/EN ISO	10993-11	2018	Valutazione biologica	Biological evaluation of medical devices – Part 11: Tests for systemic toxicity
DS/EN ISO	10993-12	2021	Valutazione biologica	Biological evaluation of medical devices – Part 12: Sample preparation and reference materials
DS/EN ISO	10993-13	2010	Valutazione biologica	Biological evaluation of medical devices – Part 13: Identification and quantification of degradation products from polymeric medicinal devices
DS/EN ISO	10993-18:2020 + A1:2023	2020/2023	Valutazione biologica	Biological evaluation of medical devices – Part 18: Chemical characterisation of materials
DS/ISO/TS	10993-19	2020	Valutazione biologica	Biological evaluation of medical devices – Part 19: Physico-chemical, morphological and topographical characterisation of materials
DS/EN/TR	10993-33	2015	Valutazione biologica	Biological evaluation of medical devices - Part 33: Guidance on tests to evaluate genotoxicity - Supplement to ISO 10993-3
DS/EN/TR	21726	2019	Valutazione biologica	Biological evaluation of medical devices - Application of the threshold of toxicological concern (TTC) for assessing biocompatibility of medical devices constituents
DS/EN	12469	2000	Camere bianche	Biotechnology – Performance criteria for microbiological safety cabinets
EN ISO	14644-1	2015	Camere bianche	Cleanrooms and associated controlled environments – Part 1: Classification of air cleanliness
EN ISO	14644-2	2015	Camere bianche	Cleanrooms and associated controlled environments – Part 2: Specification for Testing and Monitoring to Prove Continued Compliance with ISO 14644-1
EN ISO	14644-3:2019 + Cor. 2020-06	2019/2020	Camere bianche	Cleanrooms and associated controlled environments – Part 3: Test methods
EN ISO	14644-4	2022	Camere bianche	Cleanrooms and associated controlled environments – Part 4: Design, Construction and Start-up
EN ISO	14644-5	2004	Camere bianche	Cleanrooms and associated controlled environments – Part 5: Drift
EN ISO	14644-7	2004	Camere bianche	Cleanrooms and associated controlled environments – Part 7: Separative devices (clean air hood, gloveboxes, isolators and mini-environments)
EN ISO	14644-8	2022	Camere bianche	Cleanrooms and associated controlled environments – Part 8: Classification of airborne molecular contamination
DS/EN	17141	2020	Camere bianche	Cleanrooms and associated controlled environments – Biocontamination control
DS/EN ISO	14155	2020	Valutazione clinica	Clinical investigation of medical devices for human subjects – Good clinical practice
EN ISO	14630	2013	Sviluppo	Non-active surgical implants – General requirements
EN ISO	14971:2019 + A11:2021	2019/2021	Sviluppo	Medical devices – Application of risk management to medical devices
DS/EN	62366-1:2015 +AC:2015 +AC:2016 +A1:2020	2015/2016/2020	Sviluppo	Medical devices – Part 1: Application of usability engineering to medical devices
DS/ISO	639	2023	Etichettatura; nomenclatura; simboli	Code for individual languages and language groups
DS/ISO IEC	7810	2019	Etichettatura; nomenclatura; simboli	Identification cards – Physical characteristics
DS/EN ISO	15223-1	2021	Etichettatura; nomenclatura; simboli	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements

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Origine	N.	Anno	Oggetto	Titolo
Farmacopea europea (Ph.Eur)	1238	2024	Fabbricazione e validazione	Air Medicinal
DS/ISO	2859-1:1999 + Amd.1:2011	1999/2011	Fabbricazione e validazione	Sampling procedures for inspection by attributes • Part 1: Sampling schemes indexed by acceptable quality limit (AQL) for lot-by-lot inspection
ISO	4180	2019	Imballaggi	Packaging - Complete, filled transport packages - General rules for the compilation of performance test schedules
DS/EN ISO	11607-1:2020 + A11:2022 + A1:2023	2020/2022/2023	Imballaggi	Packaging for terminally sterilised medical devices – Part 1: Requirements for materials, sterile barrier systems and packaging systems
EN ISO	11607-2:2020 + A11:2022 + A1:2023	2020/2022/2023	Imballaggi	Packaging for terminally sterilised medical devices – Part 2: Validation requirements for forming, sealing and assembly processes
ANSI/AAMI	ST67	2019	Sterilizzazione	Sterilization of health care products - Requirements and guidance for selecting a sterility assurance level (SAL) for products labeled “sterile”
ASTM	F 1929-15	2015	Imballaggi	Standard Test Method for Detecting Seal Leaks in Porous Medical Packaging by Dye Penetration
ASTM	F 1980-21	2021	Imballaggi	Standard Guide for accelerated Aging of Sterile Barrier Systems and Medical Devices
ASTM	F 2096 – 11 (2019)	2019	Imballaggi	Standard Test Method for Detecting Gross Leaks in Medical Packaging by Internal Pressurisation (Bubble Test)
ASTM	F 1886/F1886M-16	2016	Imballaggi	Standard Test Method for Determining Integrity of Seals for Flexible Packaging by Visual Inspection
ASTM	F88/F 88M-23	2023	Imballaggi	Standard Test Method for Seal Strength of Flexible Barrier Materials
DS/EN ISO	13485:2016 + AC:2018 + A11:2021	2016/2018/2021	Garanzia di qualità (QA)-generale	Medical devices – Quality management systems – Requirements for regulatory purposes
ISO	28594	2017	Procedure di campionamento	Combined accept-zero sampling systems and process control procedures for product acceptance
ISO	2859-1:1999 + Cor.1:2001 + Amd 1:2011	1999/2001/2011	Procedure di campionamento	Sampling procedures for inspection by attributes – Part 1: Sampling schemes indexed by acceptable quality limit (AQL) for lot-by-lot inspection
DS/EN/ ISO	11138-1	2017	Sterilizzazione	Sterilisation of health care products - biological indicators – Part 1
DS/EN ISO	11138-3	2017	Sterilizzazione	Sterilisation of health care products - biological indicators – Part 3: Biological indicators for moist heat sterilisation indicators
DS/EN ISO	11138-7	2019	Sterilizzazione	Sterilization of health care products - biological indicators – Part 7: Guidance for selection, use and interpretation of results
EN ISO	17665-1	2006	Sterilizzazione	Sterilisation of health care products - Moist heat – Part 1: Requirements for the development, validation and routine control of a sterilisation process for medical devices
EN	556-1:2001 + AC:2006	2001/2006	Sterilizzazione	Sterilisation of medical devices to be designated “STERILE” – Part 1: Requirements for terminally sterilised medical devices

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Abbreviations	6
Summary of safety and clinical performance - Part I - Professional users	7
1. Device identification and general information	7
1.1 Device trade name	7
1.2 Manufacturer	7
1.3 Manufacturer single registration number (SRN)	7
1.4 Basic UDI-DI	7
1.5 Medical device nomenclature description	7
1.6 Class of device	8
1.7 Year when the first CE certificate was issued covering the device	8
1.8 Notified body name and identification number	8
2. Intended purpose, indication, target population and contraindications	8
2.1 Intended purpose, indication and target population	8
2.2 Contraindications and/or limitations	8
3. Device description	9
3.1 Description of device	9
3.1.1 Mode of action	9
3.2 Previous generations of the device	9
3.3 Any other devices or products intended to be used in combination with the device	9
4. Residual risks and warnings	9
4.1 Residual risks and undesirable effects	9
4.2 Warnings and precautions	10
4.2.1 Warnings	10
4.2.2 Precautions	10
4.3 Other information	11
5. Summary of clinical evaluation and post-market clinical follow-up	11
5.1 Equivalence	11
5.2 Clinical investigations prior to CE mark	11
5.2.1 "Proof of concept"-study	12
5.2.1.1 Clinical investigation design	12
5.2.1.2 Results	14
5.2.1.3 Conclusion	15

5.2.2	Safety assessment study	16
5.2.2.1	Clinical investigation design	16
5.2.2.2	Results	18
5.2.2.3	Conclusion	19
5.2.3	Histopathological report	20
5.2.3.1	Clinical investigation design	20
5.2.3.2	Results	20
5.2.3.3	Conclusion	23
5.2.4	Prospective, open-label, clinical investigation of intra-articular polyacrylamide hydrogel injection in subjects with knee OA ("IDA") (CIV ID: CIV-19-07-029146)	23
5.2.4.1	Clinical investigation design	24
5.2.4.2	Results	26
5.2.4.3	Conclusion	35
5.2.5	A multi-centre randomised, controlled, double-blind clinical investigation of intra-articular polyacrylamide hydrogel injection in subjects with knee osteoarthritis ("ROSA") (CIV ID: CIV-19-01-026964)	35
5.2.5.1	Clinical investigation design	36
5.2.5.2	Results	38
5.2.5.3	Conclusion	47
5.3	Other clinical data	48
5.4	Summary of the clinical performance and safety	48
5.5	Post-market clinical follow-up	49
6.	Possible therapeutic alternatives	50
6.1	Intra-articular treatments	50
6.1.1	Corticosteroids	50
6.1.2	Hyaluronic acid	50
6.1.3	Noltrex	51
6.2	Surgery	51
7.	Suggested profile and training for users	52
8.	Harmonised standards and applied Common Specifications	52
9.	Revision history	55
10.	Bibliography	56
	Samenvatting van veiligheid en klinische prestaties - Deel II – Patiënten	62
1.	Hulpmiddelidentificatie en algemene gegevens	62
1.1	Handelsnaam van het hulpmiddel	62
1.2	Fabrikant	62

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Table 8. Patient characteristics.

Patient	Age / Sex	Disease	OA grade	KL score*	Previous treatment	“Hydrogel B” volume	Time since “Hydrogel B” injection (months)
Case 1	35 years - female	Injury	1	-	Surgery	5 + 3 mL	18
Case 2	84 years - female	OA	2	3	Steroid	3 + 3 mL	5
Case 3	77 years - female	OA	2	-	None	5 + 5 mL	29
Case 4	69 years - female	OA	2	2	None	3 mL	33
Case 5	67 years - female	OA	3	3	None	3 + 3 mL	7
Case 6	74 years - female	OA	3	2	None	3 + 3 mL	5
Case 7	82 years - male	OA	2	3	Steroid	3 + 3 mL	16

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In essence, a similar histological pattern was seen in all seven cases. “Hydrogel B” was found to be integrated in the synovial membrane, and outer synovial lining cells had entered the gel and established a de novo lining layer. The integrated gel appears to contribute with a considerable augmentation effect to the inner side of the capsule, facing the cavity. Thus, extensive papillary folding of the synovial membrane was seen in all cases except one (case 1). In one knee (case 1) the

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week 26, week 52, year 2, year 3 and year 4 are summarised for the ITT analysis set together with the respective changes from baseline in Table 13. The changes from baseline are shown graphically in Figure 6 for the ITT analysis set, both as observed and using the BOCF approach, and for the year-2 extension subset.

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Descriptive results:

The scores for the transformed WOMAC pain subscale at baseline, weeks 4, 12, 26 and 52 are summarised for the FAS analysis set together with the respective changes from baseline in [Table 14](#). The mean transformed WOMAC pain score was lowest for both arms at 12 weeks. At 26 and 52 weeks the mean transformed WOMAC pain score appeared to increase slightly in the Synvisc-one arm compared to a more stable response in the Arthroamid® arm.

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The primary analysis (FAS) showed a treatment difference at 26 weeks of 3.59 (95% CI: -0.87; 8.05) in favour of the Arthrosamid® arm over the Synvisc-one arm (Table 15). The non-inferiority objective was met, but superiority was not reached.

For the ITT and PP analysis sets, the treatment differences were also in favour of the Arthrosamid® arm over the Synvisc-one arm (Table 15) and the non-inferiority objective was met.

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WOMAC pain subscale

The analysis of the change from baseline in the transformed WOMAC pain subscale (0-100) from baseline to weeks 4, 12 and 52 was performed using the same MMRM as for the primary endpoint. The analysis showed no significant difference in the change from baseline transformed WOMAC pain subscale scores between the Arthrosamid® and Synvisc-one arms (Table 16), although the difference at week 52 approached significance (p=0.0572).

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WOMAC stiffness subscale

The analysis of the change from baseline in the transformed WOMAC stiffness subscale (0-100) was performed using the same MMRM as for the primary endpoint. The analysis showed a statistically significant treatment difference at week 4 in favour of Synvisc-one (p = 0.0315) while treatment differences in favour of Arthroamid® at weeks 12, 26 and 52 were not significant ([Table 17](#)).

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WOMAC physical function subscale

The analysis of the change from baseline in the transformed WOMAC physical function subscale (0-100) was performed using the same MMRM as for the primary endpoint. The analysis showed a treatment difference in favour of Synvisc-one at week 4 and treatment differences in favour of Arthroamid® at weeks 12, 26 and 52. None of the differences were statistically significant ([Table 18](#)).

Table 18. Analysis of change from baseline in transformed WOMAC physical function subscale (0-100) - ITT analysis set.

	N	LSMean (95% CI)	Treatment difference (95% CI)	p-value
Week 4				
Synvisc-One	120	-17.22 (-20.29; -14.15)		
Arthroamid®	119	-13.46 (-16.54; -10.37)	-3.77 (-8.12; 0.59)	0.0896
Week 12				
Synvisc-One	119	-18.12 (-20.91; -15.32)		
Arthroamid®	115	-21.36 (-24.19; -18.54)	3.25 (-0.73; 7.22)	0.1091
Week 26				
Synvisc-One	117	-18.13 (-21.30; -14.96)		
Arthroamid®	114	-18.90 (-22.10; -15.71)	0.77 (-3.73; 5.27)	0.7357
Week 52				
Synvisc-One	109	-15.21 (-18.60; -11.82)		
Arthroamid®	107	-17.74 (-21.16; -14.33)	2.54 (-2.28; 7.35)	0.3006

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The analysis of the change from baseline in the patient global assessment was performed using the same MMRM as for the primary endpoint. The analysis showed a treatment difference in favour of Synvisc-one at week 4 and treatment differences in favour of Arthroamid® at weeks 12, 26 and 52. None of the differences were statistically significant ([Table 19](#)).

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Noltrex is based on cross-linked polyacrylamide (4.0% ± 1.5%) with silver ions added (0.0001 0.0025%). Intra-articularly injected Noltrex is claimed to cover the entire joint surface and the synovial lining with a uniform layer of lubricating gel, creating a shock-absorbing effect which helps to preserve the integrity of cartilage/stop further destruction of the cartilage. The joint space is claimed to widen due to viscosupplementation, leading to mechanical protection of articular cartilages. The cross-linked polymer is claimed to have a long-lasting effect compared to HA-based products. However, as opposed to Arthrosamid®, Noltrex is eliminated from the body, and the injections are to be repeated every 6 months - 1.5 years.

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Origin	No.	Year	Subject	Title
DS/EN ISO	10993-18:2020 + A1:2023	2020/2023	Biological evaluation	Biological evaluation of medical devices – Part 18: Chemical characterisation of materials
DS/ISO/TS	10993-19	2020	Biological evaluation	Biological evaluation of medical devices – Part 19: Physico-chemical, morphological and topographical characterisation of materials
DS/EN/TR	10993-33	2015	Biological evaluation	Biological evaluation of medical devices - Part 33: Guidance on tests to evaluate genotoxicity - Supplement to ISO 10993-3
DS/EN/TR	21726	2019	Biological evaluation	Biological evaluation of medical devices - Application of the threshold of toxicological concern (TTC) for assessing biocompatibility of medical devices constituents
DS/EN	12469	2000	Cleanrooms	Biotechnology - Performance criteria for microbiological safety cabinets
EN ISO	14644-1	2015	Cleanrooms	Cleanrooms and associated controlled environments – Part 1: Classification of air cleanliness
EN ISO	14644-2	2015	Cleanrooms	Cleanrooms and Associated Controlled Environments - Part 2: Specification for Testing and Monitoring to Prove Continued Compliance with ISO 14644-1
EN ISO	14644-3:2019 + Cor. 2020-06	2019/2020	Cleanrooms	Cleanrooms and Associated Controlled Environments - Part 3: Test methods
EN ISO	14644-4	2022	Cleanrooms	Cleanrooms and Associated Controlled Environments - Part 4: Design, Construction and Start-up
EN ISO	14644-5	2004	Cleanrooms	Cleanrooms and Associated Controlled Environments - Part 5: Drift
EN ISO	14644-7	2004	Cleanrooms	Cleanrooms and Associated Controlled Environments - Part 7: Separative devices (clean air hood, gloveboxes, isolators and mini-environments)
EN ISO	14644-8	2022	Cleanrooms	Cleanrooms and Associated Controlled Environments - Part 8: Classification of airborne molecular contamination
DS/EN	17141	2020	Cleanrooms	Cleanrooms and associated controlled environments - Bio-contamination control
DS/EN ISO	14155	2020	Clinical evaluation	Clinical investigation of medical devices for human subjects – Good clinical practice
EN ISO	14630	2013	Development	Non-active surgical implants – General requirements
EN ISO	14971:2019 + A11:2021	2019/2021	Development	Medical devices – Application of risk management to medical devices
DS/EN	62366-1:2015 +AC:2015 +AC:2016 +A1:2020	2015/2016/2020	Development	Medical Devices - Part 1: Application of usability engineering to medical devices
DS/ISO	639	2023	Labelling; Nomenclature; Symbols	Code for individual languages and language groups
DS/ISO IEC	7810	2019	Labelling; Nomenclature; Symbols	Identification cards - Physical characteristics
DS/EN ISO	15223-1	2021	Labelling; Nomenclature; Symbols	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements
DS/EN ISO	15223-1 Bil.1	2021	Labelling; Nomenclature; Symbols	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements - Fact Sheet 1: Graphical symbols in electronic format
DS/ISO	15223-2	2011	Labelling; Nomenclature; Symbols	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied - Part 2: Symbol Development, selection and validation
DS/EN ISO	20417	2021	Labelling; Nomenclature; Symbols	Information supplied by the manufacturer with medical devices
USP	61	2023	Manufacturing and validation	Microbiological Examination of Non-Sterile Products: microbial enumeration tests labs and test for specified organisms
USP	85	2023	Manufacturing and validation	Bacterial Endotoxins Test
USP	87	2023	Manufacturing and validation	Biological Reactivity Tests, In vitro

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23. Sun, S.F., Hsu, C.W., Lin, H.S., Liou, I.H., Chen, Y.H., Hung, C.L., 2017. Comparison of Single Intra-Articular Injection of Novel Hyaluronan (HYA-JOINT Plus) with Synvisc-One for Knee Osteoarthritis: A Randomized, Controlled, Double-Blind Trial of Efficacy and Safety. *The Journal of bone and joint surgery. American volume* 99, 462-471.
24. Zhao, H., Liu, H., Liang, X., Li, Y., Wang, J., Liu, C., 2016. Hylan G-F 20 Versus Low Molecular Weight Hyaluronic Acids for Knee Osteoarthritis: A Meta-Analysis. *BioDrugs : clinical immunotherapeutics, biopharmaceuticals and gene therapy* 30, 387-396.
25. Balazs, E.A., 2004. Viscosupplementation for treatment of osteoarthritis: from initial discovery to current status and results. *Surgical technology international* 12, 278-289.
26. He, W.W., Kuang, M.J., Zhao, J., Sun, L., Lu, B., Wang, Y., Ma, J.X., Ma, X.L., 2017. Efficacy and safety of intraarticular hyaluronic acid and corticosteroid for knee osteoarthritis: A meta-analysis. *International journal of surgery (London, England)* 39, 95-103.
27. Leighton, R., Akemark, C., Therrien, R., Richardson, J.B., Andersson, M., Todman, M.G., Arden, N.K., 2014. NASHA hyaluronic acid vs. methylprednisolone for knee osteoarthritis: a prospective, multi-centre, randomized, non-inferiority trial. *Osteoarthritis and cartilage / OARS, Osteoarthritis Research Society* 22, 17-25.
28. Phillips, M., Vannabouathong, C., Devji, T., Patel, R., Gomes, Z., Patel, A., Dixon, M., Bhandari, M., 2020. Differentiating factors of intra-articular injectables have a meaningful impact on knee osteoarthritis outcomes: a network meta-analysis. *Knee surgery, sports traumatology, arthroscopy: official journal of the ESSKA* 28, 3031-3039.
29. Leighton, R., Fitzpatrick, J., Smith, H., Crandall, D., Flannery, C.R., Conrozier, T., 2018. Systematic clinical evidence review of NASHA (Durolane hyaluronic acid) for the treatment of knee osteoarthritis. *Open Access Rheumatol* 10, 43-54.
30. Maheu, E., Bannuru, R.R., Herrero-Beaumont, G., Allali, F., Bard, H., Migliore, A., 2019. Why we should definitely include intra-articular hyaluronic acid as a therapeutic option in the management of knee osteoarthritis: Results of an extensive critical literature review. *Seminars in arthritis and rheumatism* 48, 563-572.
31. Maheu, E., Rannou, F., Reginster, J.Y., 2016. Efficacy and safety of hyaluronic acid in the management of osteoarthritis: Evidence from real-life setting trials and surveys. *Seminars in arthritis and rheumatism* 45, S28-33.
32. Tiendrebeogo, E., Choueiri, M., Chevalier, X., Conrozier, T., Eymard, F., 2021. Does the Presence of Neuropathic Pain Influence the Response to Hyaluronic Acid in Patients with Knee Osteoarthritis? *Cartilage* 13, 1548S-1556S.
33. Singh, H., Knapik, D.M., Polce, E.M., Eikani, C.K., Bjornstad, A.H., Gursoy, S., Perry, A.K., Westrick, J.C., Yanke, A.B., Verma, N.N., Cole, B.J., Chahla, J.A., 2022. Relative Efficacy of Intra-articular Injections in the Treatment of Knee Osteoarthritis: A Systematic Review and Network Meta-analysis. *Am J Sports Med* 50, 3140-3148.

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Opphav	Nr.	År	Tema	Tittel
DS/EN ISO	10993-1	2020	Biologisk evaluering	Biological evaluation of medical devices – Part 1 Evaluation and testing
DS/EN ISO	10993-3	2014	Biologisk evaluering	Biological evaluation of medical devices – Part 3: Test for genotoxicity, carcinogenicity and reproductive toxicity
DS/EN ISO	10993-5	2009	Biologisk evaluering	Biological evaluation of medical devices – Part 5: Tests for in vitro cytotoxicity
DS/EN ISO	10993-6	2016	Biologisk evaluering	Biological evaluation of medical devices – Part 6: Tests for local effects after implantation
DS/EN ISO	10993-9	2021	Biologisk evaluering	Biological evaluation of medical devices – Part 9 Framework for identification and degradation of potential degradation products
DS/EN ISO	10993-10	2023	Biologisk evaluering	Biological evaluation of medical devices – Part 10: Tests for irritation and delayed-type hypersensitivity
DS/EN ISO	10993-11	2018	Biologisk evaluering	Biological evaluation of medical devices – Part 11: Tests for systemic toxicity
DS/EN ISO	10993-12	2021	Biologisk evaluering	Biological evaluation of medical devices - Part 12: Sample preparation and reference materials
DS/EN ISO	10993-13	2010	Biologisk evaluering	Biological evaluation of medical devices – Part 13: Identification and quantification of degradation products from polymeric medicinal devices
DS/EN ISO	10993-18:2020 + A1:2023	2020/2023	Biologisk evaluering	Biological evaluation of medical devices – Part 18: Chemical characterisation of materials
DS/ISO/TS	10993-19	2020	Biologisk evaluering	Biological evaluation of medical devices – Part 19: Physico-chemical, morphological and topographical characterisation of materials
DS/EN/TR	10993-33	2015	Biologisk evaluering	Biological evaluation of medical devices – Part 33: Guidance on tests to evaluate genotoxicity - Supplement to ISO 10993-3
DS/EN/TR	21726	2019	Biologisk evaluering	Biological evaluation of medical devices – Application of the threshold of toxicological concern (TTC) for assessing biocompatibility of medical devices constituents
DS/EN	12469	2000	Rentrom	Biotechnology – Performance criteria for microbiological safety cabinets
EN ISO	14644-1	2015	Rentrom	Cleanrooms and associated controlled environments – Part 1: Classification of air cleanliness
EN ISO	14644-2	2015	Rentrom	Cleanrooms and Associated Controlled Environments – Part 2: Specification for Testing and Monitoring to Prove Continued Compliance with ISO 14644-1
EN ISO	14644-3:2019 + Cor. 2020-06	2019/2020	Rentrom	Cleanrooms and Associated Controlled Environments – Part 3: Test methods
EN ISO	14644-4	2022	Rentrom	Cleanrooms and Associated Controlled Environments – Part 4: Design, Construction and Start-up
EN ISO	14644-5	2004	Rentrom	Cleanrooms and Associated Controlled Environments – Part 5: Drift
EN ISO	14644-7	2004	Rentrom	Cleanrooms and Associated Controlled Environments - Part 7: Separative devices (clean air hood, gloveboxes, isolators and mini-environments)
EN ISO	14644-8	2022	Rentrom	Cleanrooms and Associated Controlled Environments – Part 8: Classification of airborne molecular contamination
DS/EN	17141	2020	Rentrom	Cleanrooms and associated controlled environments – Biocontamination control
DS/EN ISO	14155	2020	Klinisk evaluering	Clinical investigation of medical devices for human subjects – Good clinical practice
EN ISO	14630	2013	Utvikling	Non-active surgical implants – General requirements
EN ISO	14971:2019 + A11:2021	2019/2021	Utvikling	Medical devices – Application of risk management to medical devices
DS/EN	62366-1:2015 +AC:2015 +AC:2016 +A1:2020	2015/2016/2020	Utvikling	Medical Devices - Part 1: Application of usability engineering to medical devices

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Opphav	Nr.	År	Tema	Tittel
DS/ISO	2859-1:1999 + Amd.1:2011	1999/2011	Produksjon og validering	Sampling procedures for inspection by attributes • Part 1: Sampling schemes indexed by acceptable quality limit (AQL) for lot-by-lot inspection
ISO	4180	2019	Emballasje	Packaging – Complete, filled transport packages – General rules for the compilation of performance test schedules
DS/EN ISO	11607-1:2020 + A11:2022 + A1:2023	2020/2022/2023	Emballasje	Packaging for terminally sterilised medical devices – Part 1: Requirements for materials, sterile barrier systems and packaging systems
EN ISO	11607-2:2020 + A11:2022 + A1:2023	2020/2022/2023	Emballasje	Packaging for terminally sterilised medical devices – Part 2: Validation requirements for forming, sealing and assembly processes
ANSI/AAMI	ST67	2019	Sterilisering	Sterilization of health care products – requirements and guidance for selecting a sterility assurance level (SAL) for products labeled "sterile"
ASTM	F 1929-15	2015	Emballasje	Standard Test Method for Detecting Seal Leaks in Porous Medical Packaging by Dye Penetration
ASTM	F 1980-21	2021	Emballasje	Standard Guide for accelerated Aging of Sterile Barrier Systems and Medical Devices
ASTM	F 2096 – 11 (2019)	2019	Emballasje	Standard Test Method for Detecting Gross Leaks in Medical Packaging by Internal Pressurisation (Bubble Test)
ASTM	F 1886/F1886M-16	2016	Emballasje	Standard Test Method for Determining Integrity of Seals for Flexible Packaging by Visual Inspection
ASTM	F88/F 88M-23	2023	Emballasje	Standard Test Method for Seal Strength of Flexible Barrier Materials
DS/EN ISO	13485:2016 + AC:2018 + A11:2021	2016/2018/2021	Generell kvalitetskontroll	Medical devices – Quality management systems – Requirements for regulatory purposes
ISO	28594	2017	Prøvetakingsprosedyrer	Combined accept-zero sampling systems and process control procedures for product acceptance
ISO	2859-1:1999 + Cor.1:2001 + Amd.1:2011	1999/2001/2011	Prøvetakingsprosedyrer	Sampling procedures for inspection by attributes - Part 1: Sampling schemes indexed by acceptable quality limit (AQL) for lot-by-lot inspection
DS/EN/ ISO	11138-1	2017	Sterilisering	Sterilisation of health care products – biological indicators – Part 1
DS/EN ISO	11138-3	2017	Sterilisering	Sterilisation of health care products – biological indicators – Part 3 Biological indicators for moist heat sterilisation indicators
DS/EN ISO	11138-7	2019	Sterilisering	Sterilization of health care products – biological indicators – Part 7: Guidance for selection, use and interpretation of results
EN ISO	17665-1	2006	Sterilisering	Sterilisation of health care products – Moist heat – Part 1: Requirements for the development, validation and routine control of a sterilisation process for medical devices
EN	556-1:2001 + AC:2006	2001/2006	Sterilisering	Sterilisation of medical devices to be designated "STERILE" – Part 1: Requirements for terminally sterilised medical devices

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